

产品名称: CD321 Mouse Monoclonal Antibody
产品货号: AMM82612



产品概述 (Summary)

产品名称 (Production Name)	CD321 Mouse Monoclonal Antibody
描述 (Description)	Mouse monoclonal Antibody
宿主 (Host)	Mouse
应用 (Application)	IHC,ELISA,FC
种属反应性 (Reactivity)	Human

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	Mouse IgG1
克隆 (Clonality)	Monoclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	Purified antibody in PBS with 0.05% sodium azide
纯化方式 (Purification)	Affinity Purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	CD321
别名 (Alternative Names)	F11R;JAM; KAT; JAM1; JAMA; JCAM; PAM-1
基因 ID (Gene ID)	50848.0
蛋白 ID (SwissProt ID)	Q9Y624.Purified recombinant fragment of human CD321 (AA: extra 28-238) expressed in E. Coli.

产品应用 (Application)

稀释比 (Dilution Ratio)	IHC 1:200-1:1000,ELISA 1:5000-1:20000,FC 1:200-1:400
蛋白分子量 (Molecular Weight)	32.6kDa

研究背景 (Background)

Tight junctions represent one mode of cell-to-cell adhesion in epithelial or endothelial cell sheets, forming continuous seals


产品名称: CD321 Mouse Monoclonal Antibody
产品货号: AMM82612




around cells and serving as a physical barrier to prevent solutes and water from passing freely through the paracellular space. The protein encoded by this immunoglobulin superfamily gene member is an important regulator of tight junction assembly in epithelia. In addition, the encoded protein can act as (1) a receptor for reovirus, (2) a ligand for the integrin LFA1, involved in leukocyte transmigration, and (3) a platelet receptor. Multiple 5' alternatively spliced variants, encoding the same protein, have been identified but their biological validity has not been established.

研究领域 (Research Area)

图片 (Image Data)




Black line: Control Antigen (100 ng);Purple line: Antigen (10ng); Blue line: Antigen (50 ng); Red line:Antigen (100 ng)




Flow cytometric analysis of HT-29 cells using CD321 mouse mAb (green) and negative control (red).

产品名称: CD321 Mouse Monoclonal Antibody
产品货号: AMM82612



Immunohistochemical analysis of paraffin-embedded human lung cancer tissues using CD321 mouse mAb with DAB staining.



Immunohistochemical analysis of paraffin-embedded human liver cancer tissues using CD321 mouse mAb with DAB staining.

注意事项 (Note)

For research use only .