

产品名称: KRT2 Mouse Monoclonal Antibody
产品货号: AMM82658



产品概述 (Summary)

产品名称 (Production Name)	KRT2 Mouse Monoclonal Antibody
描述 (Description)	Mouse monoclonal Antibody
宿主 (Host)	Mouse
应用 (Application)	WB,IHC,ICC,ELISA,FC
种属反应性 (Reactivity)	Human

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	Mouse IgG2a
克隆 (Clonality)	Monoclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	Purified antibody in PBS with 0.05% sodium azide
纯化方式 (Purification)	Affinity Purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	KRT2
别名 (Alternative Names)	K2e; KRTE; CK-2e; KRT2A; KRT2E
基因 ID (Gene ID)	3849.0
蛋白 ID (SwissProt ID)	P35908.Purified recombinant fragment of human KRT2 (AA: 177-490) expressed in E. Coli.

产品应用 (Application)

稀释比 (Dilution Ratio)	WB 1:500-1:2000,IHC 1:200-1:1000,ICC 1:200-1:1000,ELISA 1:5000-1:20000,FC 1:200-1:400
蛋白分子量 (Molecular Weight)	65.4kDa

研究背景 (Background)


产品名称: KRT2 Mouse Monoclonal Antibody
产品货号: AMM82658




The protein encoded by this gene is a member of the keratin gene family. The type II cytokeratins consist of basic or neutral proteins which are arranged in pairs of heterotypic keratin chains coexpressed during differentiation of simple and stratified epithelial tissues. This type II cytokeratin is expressed largely in the upper spinous layer of epidermal keratinocytes and mutations in this gene have been associated with bullous congenital ichthyosiform erythroderma. The type II cytokeratins are clustered in a region of chromosome 12q12-q13.

研究领域 (Research Area)

图片 (Image Data)




Black line: Control Antigen (100 ng);Purple line: Antigen (10ng); Blue line: Antigen (50 ng); Red line:Antigen (100 ng)




Western blot analysis using KRT2 mouse mAb against A431 (1) cell lysate.


产品名称: KRT2 Mouse Monoclonal Antibody
产品货号: AMM82658



Immunofluorescence analysis of HeLa cells using KRT2 mouse mAb (green). Blue: DRAQ5 fluorescent DNA dye. Red: Actin filaments have been labeled with Alexa Fluor- 555 phalloidin. Secondary antibody from Fisher (Cat#: 35503)




Flow cytometric analysis of HeLa cells using KRT2 mouse mAb (green) and negative control (red).



Immunohistochemical analysis of paraffin-embedded human lung cancer tissues using KRT2 mouse mAb with DAB staining.

产品名称: KRT2 Mouse Monoclonal Antibody
产品货号: AMM82658



Immunohistochemical analysis of paraffin-embedded human bladder cancer tissues using KRT2 mouse mAb with DAB staining.

注意事项 (Note)

For research use only .