

产品名称: ATP1A1 Mouse Monoclonal Antibody
产品货号: AMM82660



产品概述 (Summary)

产品名称 (Production Name)	ATP1A1 Mouse Monoclonal Antibody
描述 (Description)	Mouse monoclonal Antibody
宿主 (Host)	Mouse
应用 (Application)	WB,IHC,ICC,ELISA,FC
种属反应性 (Reactivity)	Human

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	Mouse IgG2a
克隆 (Clonality)	Monoclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	Purified antibody in PBS with 0.05% sodium azide
纯化方式 (Purification)	Affinity Purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	ATP1A1
别名 (Alternative Names)	CMT2DD; HOMGSMR2
基因 ID (Gene ID)	476.0
蛋白 ID (SwissProt ID)	P05023.Purified recombinant fragment of human ATP1A1 (AA: 153-288) expressed in E. Coli.

产品应用 (Application)

稀释比 (Dilution Ratio)	WB 1:500-1:2000,IHC 1:200-1:1000,ICC 1:200-1:1000,ELISA 1:5000-1:20000,FC 1:200-1:400
蛋白分子量 (Molecular Weight)	112.8kDa

研究背景 (Background)


产品名称: ATP1A1 Mouse Monoclonal Antibody
产品货号: AMM82660




The protein encoded by this gene belongs to the family of P-type cation transport ATPases, and to the subfamily of Na⁺/K⁺ -ATPases. Na⁺/K⁺ -ATPase is an integral membrane protein responsible for establishing and maintaining the electrochemical gradients of Na and K ions across the plasma membrane. These gradients are essential for osmoregulation, for sodium-coupled transport of a variety of organic and inorganic molecules, and for electrical excitability of nerve and muscle. This enzyme is composed of two subunits, a large catalytic subunit (alpha) and a smaller glycoprotein subunit (beta). The catalytic subunit of Na⁺/K⁺ -ATPase is encoded by multiple genes. This gene encodes an alpha 1 subunit. Multiple transcript variants encoding different isoforms have been found for this gene.

研究领域 (Research Area)

图片 (Image Data)




Black line: Control Antigen (100 ng);Purple line: Antigen (10ng); Blue line: Antigen (50 ng); Red line:Antigen (100 ng)




Western blot analysis using ATP1A1 mouse mAb against Hela (1), HepG2 (2), and A431 (3) cell lysate.


产品名称: ATP1A1 Mouse Monoclonal Antibody
产品货号: AMM82660



Immunofluorescence analysis of HeLa cells using ATP1A1 mouse mAb (green). Blue: DRAQ5 fluorescent DNA dye. Red: Actin filaments have been labeled with Alexa Fluor- 555 phalloidin. Secondary antibody from Fisher (Cat#: 35503)




Flow cytometric analysis of HeLa cells using ATP1A1 mouse mAb (green) and negative control (red).



Immunohistochemical analysis of paraffin-embedded human lung cancer tissues using ATP1A1 mouse mAb with DAB staining.

产品名称: ATP1A1 Mouse Monoclonal Antibody
产品货号: AMM82660



Immunohistochemical analysis of paraffin-embedded human rectum cancer tissues using ATP1A1 mouse mAb with DAB staining.

注意事项 (Note)

For research use only .