

产品名称: GZMK Mouse Monoclonal Antibody
产品货号: AMM82681



产品概述 (Summary)

产品名称 (Production Name)	GZMK Mouse Monoclonal Antibody
描述 (Description)	Mouse monoclonal Antibody
宿主 (Host)	Mouse
应用 (Application)	IHC,ELISA,FC
种属反应性 (Reactivity)	Human

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	Mouse IgG2b
克隆 (Clonality)	Monoclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	Purified antibody in PBS with 0.05% sodium azide
纯化方式 (Purification)	Affinity Purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	GZMK
别名 (Alternative Names)	TRYP2
基因 ID (Gene ID)	3003.0
蛋白 ID (SwissProt ID)	P49863.Purified recombinant fragment of human GZMK (AA: 27-264) expressed in E. Coli.

产品应用 (Application)

稀释比 (Dilution Ratio)	IHC 1:200-1:1000,ELISA 1:5000-1:20000,FC 1:200-1:400
蛋白分子量 (Molecular Weight)	29kda

研究背景 (Background)

This gene product is a member of a group of related serine proteases from the cytoplasmic granules of cytotoxic


产品名称: GZMK Mouse Monoclonal Antibody
产品货号: AMM82681




lymphocytes. Cytolytic T lymphocytes (CTL) and natural killer (NK) cells share the remarkable ability to recognize, bind, and lyse specific target cells. They are thought to protect their host by lysing cells bearing on their surface 'nonself' antigens, usually peptides or proteins resulting from infection by intracellular pathogens. The protein described here lacks consensus sequences for N-glycosylation present in other granzymes. [provided by RefSeq, Jul 2008]

研究领域 (Research Area)


图片 (Image Data)



Black line: Control Antigen (100 ng);Purple line: Antigen (10ng); Blue line: Antigen (50 ng); Red line:Antigen (100 ng)




Flow cytometric analysis of Jurkat cells using GZMK mouse mAb (green) and negative control (red).



产品名称: GZMK Mouse Monoclonal Antibody
产品货号: AMM82681



Immunohistochemical analysis of paraffin-embedded human endometrial cancer tissues using GZMK mouse mAb with DAB staining.



Immunohistochemical analysis of paraffin-embedded human tonsil tissues using GZMK mouse mAb with DAB staining.

注意事项 (Note)

For research use only .