

产品名称: CDK4 Mouse Monoclonal Antibody
产品货号: AMM82684



产品概述 (Summary)

产品名称 (Production Name)	CDK4 Mouse Monoclonal Antibody
描述 (Description)	Mouse monoclonal Antibody
宿主 (Host)	Mouse
应用 (Application)	WB,IHC,ELISA,FC
种属反应性 (Reactivity)	Human, Mouse, Monkey

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	Mouse IgG2a
克隆 (Clonality)	Monoclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	Purified antibody in PBS with 0.05% sodium azide
纯化方式 (Purification)	Affinity Purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	CDK4
别名 (Alternative Names)	CMM3; PSK-J3
基因 ID (Gene ID)	1019.0
蛋白 ID (SwissProt ID)	P11802.Purified recombinant fragment of human CDK4 (AA: 77-303) expressed in E. Coli.

产品应用 (Application)

稀释比 (Dilution Ratio)	WB 1:500-1:2000,IHC 1:200-1:1000,ELISA 1:5000-1:20000,FC 1:200-1:400
蛋白分子量 (Molecular Weight)	33.7kda

研究背景 (Background)

The protein encoded by this gene is a member of the Ser/Thr protein kinase family. This protein is highly similar to the

产品名称: CDK4 Mouse Monoclonal Antibody
产品货号: AMM82684




gene products of *S. cerevisiae* cdc28 and *S. pombe* cdc2. It is a catalytic subunit of the protein kinase complex that is important for cell cycle G1 phase progression. The activity of this kinase is restricted to the G1-S phase, which is controlled by the regulatory subunits D-type cyclins and CDK inhibitor p16(INK4a). This kinase was shown to be responsible for the phosphorylation of retinoblastoma gene product (Rb). Mutations in this gene as well as in its related proteins including D-type cyclins, p16(INK4a) and Rb were all found to be associated with tumorigenesis of a variety of cancers. Multiple polyadenylation sites of this gene have been reported.


研究领域 (Research Area)

PI3K-Akt signaling pathway

图片 (Image Data)




Black line: Control Antigen (100 ng);Purple line: Antigen (10ng); Blue line: Antigen (50 ng); Red line:Antigen (100 ng)




Western blot analysis using CDK4 mouse mAb against COS-7 (1), and NIH/3T3 (2) cell lysate.


产品名称: CDK4 Mouse Monoclonal Antibody
产品货号: AMM82684




Flow cytometric analysis of HeLa cells using CDK4 mouse mAb (green) and negative control (red).



Flow cytometric analysis of HepG2 cells using CDK4 mouse mAb (green) and negative control (red).




Flow cytometric analysis of Raji cells using CDK4 mouse mAb (green) and negative control (red).




产品名称: CDK4 Mouse Monoclonal Antibody
产品货号: AMM82684



Flow cytometric analysis of THP-1 cells using CDK4 mouse mAb (green) and negative control (red).



Immunohistochemical analysis of paraffin-embedded human cervical cancer tissues using CDK4 mouse mAb with DAB staining.



Immunohistochemical analysis of paraffin-embedded human kidney cancer tissues using CDK4 mouse mAb with DAB staining.

注意事项 (Note)

For research use only .