

产品名称: CD109 Mouse Monoclonal Antibody
产品货号: AMM82693



产品概述 (Summary)

产品名称 (Production Name)	CD109 Mouse Monoclonal Antibody
描述 (Description)	Mouse monoclonal Antibody
宿主 (Host)	Mouse
应用 (Application)	IHC, ICC, ELISA, FC
种属反应性 (Reactivity)	Human

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	Mouse IgG1
克隆 (Clonality)	Monoclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	Purified antibody in PBS with 0.05% sodium azide
纯化方式 (Purification)	Affinity Purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	CD109
别名 (Alternative Names)	p180; r150; CPAMD7
基因 ID (Gene ID)	135228.0
蛋白 ID (SwissProt ID)	Q6YHK3. Purified recombinant fragment of human CD109 (AA: extra 1274-1421) expressed in E. Coli.

产品应用 (Application)

稀释比 (Dilution Ratio)	IHC 1:200-1:1000, ICC 1:20-1:100, ELISA 1:5000-1:20000, FC 1:200-1:400
蛋白分子量 (Molecular Weight)	161.7kDa

研究背景 (Background)

This gene encodes a glycosyl phosphatidylinositol (GPI)-linked glycoprotein that localizes to the surface of platelets,


产品名称: CD109 Mouse Monoclonal Antibody
产品货号: AMM82693




activated T-cells, and endothelial cells. The protein binds to and negatively regulates signalling by transforming growth factor beta (TGF-beta). Multiple transcript variants encoding different isoforms have been found for this gene.

研究领域 (Research Area)


图片 (Image Data)



Black line: Control Antigen (100 ng);Purple line: Antigen (10ng); Blue line: Antigen (50 ng); Red line:Antigen (100 ng)




Immunofluorescence analysis of HeLa cells using CD109 mouse mAb (green). Blue: DRAQ5 fluorescent DNA dye. Red: Actin filaments have been labeled with Alexa Fluor- 555 phalloidin. Secondary antibody from Fisher (Cat#: 35503)




Flow cytometric analysis of Jurkat cells using CD109 mouse mAb (green) and negative control (red).

产品名称: CD109 Mouse Monoclonal Antibody
产品货号: AMM82693



Immunohistochemical analysis of paraffin-embedded human lung tissues using CD109 mouse mAb with DAB staining.



Immunohistochemical analysis of paraffin-embedded human stomach tissues using CD109 mouse mAb with DAB staining.

注意事项 (Note)

For research use only .