

产品名称: CD34 Mouse Monoclonal Antibody
产品货号: AMM82695



产品概述 (Summary)

产品名称 (Production Name)	CD34 Mouse Monoclonal Antibody
描述 (Description)	Mouse monoclonal Antibody
宿主 (Host)	Mouse
应用 (Application)	ICC,ELISA,FC
种属反应性 (Reactivity)	Human

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	Mouse IgG1
克隆 (Clonality)	Monoclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	Purified antibody in PBS with 0.05% sodium azide
纯化方式 (Purification)	Affinity Purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	CD34
别名 (Alternative Names)	CD34
基因 ID (Gene ID)	947.0
蛋白 ID (SwissProt ID)	P28906.Purified recombinant fragment of human CD34 (AA: extra 32-290) expressed in E. Coli.

产品应用 (Application)

稀释比 (Dilution Ratio)	ICC 1:20-1:100,ELISA 1:5000-1:20000,FC 1:200-1:400
蛋白分子量 (Molecular Weight)	40.7kDa

研究背景 (Background)

The protein encoded by this gene may play a role in the attachment of stem cells to the bone marrow extracellular matrix


产品名称: CD34 Mouse Monoclonal Antibody
产品货号: AMM82695




or to stromal cells. This single-pass membrane protein is highly glycosylated and phosphorylated by protein kinase C. Two transcript variants encoding different isoforms have been found for this gene.

研究领域 (Research Area)


图片 (Image Data)



Black line: Control Antigen (100 ng);Purple line: Antigen (10ng); Blue line: Antigen (50 ng); Red line:Antigen (100 ng)



Immunofluorescence analysis of HeLa cells using CD34 mouse mAb (green). Blue: DRAQ5 fluorescent DNA dye. Red: Actin filaments have been labeled with Alexa Fluor- 555 phalloidin. Secondary antibody from Fisher (Cat#: 35503)



Flow cytometric analysis of Jurkat cells using CD34 mouse mAb (green) and negative control (red).

产品名称: CD34 Mouse Monoclonal Antibody
产品货号: AMM82695



注意事项 (Note)

For research use only .