

产品名称: CXCL9 Mouse Monoclonal Antibody
产品货号: AMM82700



产品概述 (Summary)

产品名称 (Production Name)	CXCL9 Mouse Monoclonal Antibody
描述 (Description)	Mouse monoclonal Antibody
宿主 (Host)	Mouse
应用 (Application)	ELISA,FC
种属反应性 (Reactivity)	Human

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	Mouse IgG2b
克隆 (Clonality)	Monoclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	Purified antibody in PBS with 0.05% sodium azide
纯化方式 (Purification)	Affinity Purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	CXCL9
别名 (Alternative Names)	CMK; MIG; Humig; SCYB9; crg-10
基因 ID (Gene ID)	4283.0
蛋白 ID (SwissProt ID)	Q07325.Purified recombinant fragment of human CXCL9 (AA: 23-125) expressed in E. Coli.

产品应用 (Application)

稀释比 (Dilution Ratio)	ELISA 1:5000-1:20000,FC 1:200-1:400
蛋白分子量 (Molecular Weight)	14kd

研究背景 (Background)

This antimicrobial gene is part of a chemokine superfamily that encodes secreted proteins involved in immunoregulatory


产品名称: CXCL9 Mouse Monoclonal Antibody
产品货号: AMM82700




and inflammatory processes. The protein encoded is thought to be involved in T cell trafficking. The encoded protein binds to C-X-C motif chemokine 3 and is a chemoattractant for lymphocytes but not for neutrophils. [provided by RefSeq, Aug 2020]

研究领域 (Research Area)


图片 (Image Data)



Black line: Control Antigen (100 ng);Purple line: Antigen (10ng); Blue line: Antigen (50 ng); Red line:Antigen (100 ng)



Flow cytometric analysis of K562 cells using CXCL9 mouse mAb (green) and negative control (red).



Flow cytometric analysis of THP-1 cells using CXCL9 mouse mAb (green) and negative control (red).

产品名称: CXCL9 Mouse Monoclonal Antibody
产品货号: AMM82700



注意事项 (Note)

For research use only .