

产品名称: TSPAN33 Mouse Monoclonal Antibody
产品货号: AMM82707



产品概述 (Summary)

产品名称 (Production Name)	TSPAN33 Mouse Monoclonal Antibody
描述 (Description)	Mouse monoclonal Antibody
宿主 (Host)	Mouse
应用 (Application)	IHC,ELISA,FC
种属反应性 (Reactivity)	Human

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	Mouse IgG2b
克隆 (Clonality)	Monoclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	Purified antibody in PBS with 0.05% sodium azide
纯化方式 (Purification)	Affinity Purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	TSPAN33
别名 (Alternative Names)	PEN; PEN.; TSPAN-33
基因 ID (Gene ID)	340348.0
蛋白 ID (SwissProt ID)	Q86UF1 .Purified recombinant fragment of human TSPAN33 (AA:Extra(46-64)+Extra(118-235) expressed in E. Coli.

产品应用 (Application)

稀释比 (Dilution Ratio)	IHC 1:200-1:1000,ELISA 1:5000-1:20000,FC 1:200-1:400
蛋白分子量 (Molecular Weight)	32kDa

研究背景 (Background)

TSPAN33 is a Protein Coding gene. Diseases associated with TSPAN33 include Occlusion Precerebral Artery and Intracranial


产品名称: TSPAN33 Mouse Monoclonal Antibody
产品货号: AMM82707




Vasospasm. Gene Ontology annotations related to this gene include enzyme binding. An important paralog of this gene is TSPAN17.

研究领域 (Research Area)


图片 (Image Data)



Black line: Control Antigen (100 ng);Purple line: Antigen (10ng); Blue line: Antigen (50 ng); Red line:Antigen (100 ng)




Flow cytometric analysis of Jurkat cells using TSPAN33 mouse mAb (green) and negative control (red).



Flow cytometric analysis of MOLT4 cells using TSPAN33 mouse mAb (green) and negative control (red).

产品名称: TSPAN33 Mouse Monoclonal Antibody
产品货号: AMM82707



Immunohistochemical analysis of paraffin-embedded human endometrial cancer tissues using TSPAN33 mouse mAb with DAB staining.

注意事项 (Note)

For research use only .