

产品名称: CD35 Mouse Monoclonal Antibody
产品货号: AMM82737



产品概述 (Summary)

产品名称 (Production Name)	CD35 Mouse Monoclonal Antibody
描述 (Description)	Mouse monoclonal Antibody
宿主 (Host)	Mouse
应用 (Application)	IHC,ELISA,FC
种属反应性 (Reactivity)	Human

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	Mouse IgG1
克隆 (Clonality)	Monoclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	Purified antibody in PBS with 0.05% sodium azide
纯化方式 (Purification)	Affinity Purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	CD35
别名 (Alternative Names)	KN; C3BR; C4BR
基因 ID (Gene ID)	1378.0
蛋白 ID (SwissProt ID)	P17927.Purified recombinant fragment of human CD35 expressed in E. Coli.

产品应用 (Application)

稀释比 (Dilution Ratio)	IHC 1:500-1:1000,ELISA 1:5000-1:20000,FC 1:200-1:400
蛋白分子量 (Molecular Weight)	161.7kDa

研究背景 (Background)

This gene encodes a glycosyl phosphatidylinositol (GPI)-linked glycoprotein that localizes to the surface of platelets, activated T-cells, and endothelial cells. The protein binds to and negatively regulates signalling by transforming growth


产品名称: CD35 Mouse Monoclonal Antibody
产品货号: AMM82737




factor beta (TGF-beta). Multiple transcript variants encoding different isoforms have been found for this gene.

研究领域 (Research Area)

图片 (Image Data)



Black line: Control Antigen (100 ng);Purple line: Antigen (10ng); Blue line: Antigen (50 ng); Red line:Antigen (100 ng)



Flow cytometric analysis of THP-1 cells using CD35 mouse mAb (green) and negative control (red).

注意事项 (Note)

For research use only .