

产品名称: BCL2 Mouse Monoclonal Antibody
产品货号: AMM82738



产品概述 (Summary)

产品名称 (Production Name)	BCL2 Mouse Monoclonal Antibody
描述 (Description)	Mouse monoclonal Antibody
宿主 (Host)	Mouse
应用 (Application)	ELISA,FC
种属反应性 (Reactivity)	Human

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	Mouse IgG1
克隆 (Clonality)	Monoclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	Purified antibody in PBS with 0.05% sodium azide
纯化方式 (Purification)	Affinity Purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	BCL2
别名 (Alternative Names)	Bcl-2; PPP1R50
基因 ID (Gene ID)	596.0
蛋白 ID (SwissProt ID)	P10415.Purified recombinant fragment of human BCL2 (AA: 1-239) expressed in E. Coli.

产品应用 (Application)

稀释比 (Dilution Ratio)	ELISA 1:5000-1:20000,FC 1:200-1:400
蛋白分子量 (Molecular Weight)	26.3kDa

研究背景 (Background)

This gene encodes an integral outer mitochondrial membrane protein that blocks the apoptotic death of some cells such as

产品名称: BCL2 Mouse Monoclonal Antibody
产品货号: AMM82738




lymphocytes. Constitutive expression of BCL2, such as in the case of translocation of BCL2 to Ig heavy chain locus, is thought to be the cause of follicular lymphoma. Alternative splicing results in multiple transcript variants.


研究领域 (Research Area)

Apoptosis, TGF-beta signaling pathway, PI3K-Akt signaling pathway

图片 (Image Data)



Black line: Control Antigen (100 ng); Purple line: Antigen (10ng); Blue line: Antigen (50 ng); Red line: Antigen (100 ng)



Flow cytometric analysis of THP-1 cells using BCL2 mouse mAb (green) and negative control (red).

注意事项 (Note)

For research use only .