

产品名称: P4HB Mouse Monoclonal Antibody
产品货号: AMM82746



产品概述 (Summary)

产品名称 (Production Name)	P4HB Mouse Monoclonal Antibody
描述 (Description)	Mouse monoclonal Antibody
宿主 (Host)	Mouse
应用 (Application)	WB,IHC,ELISA,FC
种属反应性 (Reactivity)	Human

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	Mouse IgG2a
克隆 (Clonality)	Monoclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	Purified antibody in PBS with 0.05% sodium azide
纯化方式 (Purification)	Affinity Purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	P4HB
别名 (Alternative Names)	DSI,GIT,PDI,PHDB,PDIA1,PO4DB,PO4HB,PROHB,CLCRP1,ERBA2L,P4Hbeta
基因 ID (Gene ID)	5034.0
蛋白 ID (SwissProt ID)	P07237.Purified recombinant fragment of human P4HB (AA: 309-508) expressed in Mammalian.

产品应用 (Application)

稀释比 (Dilution Ratio)	WB 1:500-1:2000,IHC 1:200-1:1000,ELISA 1:5000-1:20000,FC 1:200-1:400
蛋白分子量 (Molecular Weight)	57.1kDa

研究背景 (Background)

This gene encodes the beta subunit of prolyl 4-hydroxylase, a highly abundant multifunctional enzyme that belongs to the


产品名称: P4HB Mouse Monoclonal Antibody
产品货号: AMM82746




protein disulfide isomerase family. When present as a tetramer consisting of two alpha and two beta subunits, this enzyme is involved in hydroxylation of prolyl residues in procollagen. This enzyme is also a disulfide isomerase containing two thioredoxin domains that catalyze the formation, breakage and rearrangement of disulfide bonds. Other known functions include its ability to act as a chaperone that inhibits aggregation of misfolded proteins in a concentration-dependent manner, its ability to bind thyroid hormone, its role in both the influx and efflux of S-nitrosothiol-bound nitric oxide, and its function as a subunit of the microsomal triglyceride transfer protein complex.

研究领域 (Research Area)

图片 (Image Data)




Black line: Control Antigen (100 ng); Purple line: Antigen (10ng); Blue line: Antigen (50 ng); Red line: Antigen (100 ng)




Western blot analysis using P4HB mouse mAb against Hela (1), PANC-1 (2), MCF-7 (3), THP-1 (4), SW620 (5), and HepG2 (6) cell lysate.


产品名称: P4HB Mouse Monoclonal Antibody
产品货号: AMM82746



Flow cytometric analysis of A375 cells using P4HB mouse mAb (green) and negative control (red).




Flow cytometric analysis of Hela cells using P4HB mouse mAb (green) and negative control (red).




Immunohistochemical analysis of paraffin-embedded human cervical cancer tissues using P4HB mouse mAb with DAB staining.

产品名称: P4HB Mouse Monoclonal Antibody
产品货号: AMM82746



Immunohistochemical analysis of paraffin-embedded human rectal cancer tissues using P4HB mouse mAb with DAB staining.



Immunohistochemical analysis of paraffin-embedded human lung cancer tissues using P4HB mouse mAb with DAB staining.

注意事项 (Note)

For research use only .