

产品名称: RABEP2 Mouse Monoclonal Antibody
产品货号: AMM82764



产品概述 (Summary)

产品名称 (Production Name)	RABEP2 Mouse Monoclonal Antibody
描述 (Description)	Mouse monoclonal Antibody
宿主 (Host)	Mouse
应用 (Application)	IHC,ELISA,FC
种属反应性 (Reactivity)	Human

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	Mouse IgG1
克隆 (Clonality)	Monoclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	Purified antibody in PBS with 0.05% sodium azide
纯化方式 (Purification)	Affinity Purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	RABEP2
别名 (Alternative Names)	FRA
基因 ID (Gene ID)	79874.0
蛋白 ID (SwissProt ID)	Q9H5N1.Purified recombinant fragment of human RABEP2 (AA: 10-200) expressed in HEK293-6e cells supernatant.

产品应用 (Application)

稀释比 (Dilution Ratio)	IHC 1:200-1:1000,ELISA 1:5000-1:20000,FC 1:200-1:400
蛋白分子量 (Molecular Weight)	63.5kDa

研究背景 (Background)

This gene encodes a Rabaptin, RAB GTPase Binding Effector Protein 2. The encoded protein Plays a role in membrane


产品名称: RABEP2 Mouse Monoclonal Antibody
产品货号: AMM82764




trafficking and in homotypic early endosome fusion. Participates in arteriogenesis by regulating vascular endothelial growth factor receptor 2/VEGFR2 cell surface expression and endosomal trafficking . By interacting with SDCCAG8, localizes to centrosomes and plays a critical role in ciliogenesis

研究领域 (Research Area)


图片 (Image Data)



Black line: Control Antigen (100 ng);Purple line: Antigen (10ng); Blue line: Antigen (50 ng); Red line:Antigen (100 ng)



Flow cytometric analysis of Hela cells using RABEP2 mouse mAb (green) and negative control (red).




Immunohistochemical analysis of paraffin-embedded human bladder cancer tissues using RABEP2 mouse mAb with DAB

产品名称: RABEP2 Mouse Monoclonal Antibody
产品货号: AMM82764



staining.



Immunohistochemical analysis of paraffin-embedded human ovarian cancer tissues using RABEP2 mouse mAb with DAB staining.

注意事项 (Note)

For research use only .