

产品名称: HSP70 Mouse Monoclonal Antibody
产品货号: AMM82775



产品概述 (Summary)

产品名称 (Production Name)	HSP70 Mouse Monoclonal Antibody
描述 (Description)	Mouse monoclonal Antibody
宿主 (Host)	Mouse
应用 (Application)	WB,IHC,ELISA,FC
种属反应性 (Reactivity)	Human, Mouse, Rat, Monkey

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	Mouse IgG2b
克隆 (Clonality)	Monoclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	Purified antibody in PBS with 0.05% sodium azide
纯化方式 (Purification)	Affinity Purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	HSP70
别名 (Alternative Names)	RY; APG-2; HSPH2; HSPA4 ; hsp70RY; HEL-S-5a; HS24/P52
基因 ID (Gene ID)	3308.0
蛋白 ID (SwissProt ID)	P34932.Purified recombinant fragment of human HSP70 (AA: 642-841) expressed in Mammalian

产品应用 (Application)

稀释比 (Dilution Ratio)	WB 1:500-1:2000,IHC 1:200-1:1000,ELISA 1:5000-1:20000,FC 1:200-1:400
蛋白分子量 (Molecular Weight)	94.3kDa


研究背景 (Background)

产品名称: HSP70 Mouse Monoclonal Antibody
产品货号: AMM82775




研究领域 (Research Area)


图片 (Image Data)



Black line: Control Antigen (100 ng);Purple line: Antigen (10ng); Blue line: Antigen (50 ng); Red line:Antigen (100 ng)




Western blot analysis using HSP70 mouse mAb against Hela (1), HepG2 (2), Hek293 (3), PC-12 (4), A549 (5), Jurkat (6), COS-7 (7), and NIH/3T3 (8) cell lysate.




Flow cytometric analysis of Hela cells using HSP70 mouse mAb (green) and negative control (red).


产品名称: HSP70 Mouse Monoclonal Antibody
产品货号: AMM82775



Flow cytometric analysis of Jurkat cells using HSP70 mouse mAb (green) and negative control (red).



Immunohistochemical analysis of paraffin-embedded human bladder Cancer tissues using HSP70 mouse mAb with DAB staining.



Immunohistochemical analysis of paraffin-embedded human rectal cancer tissues using HSP70 mouse mAb with DAB staining.

注意事项 (Note)

For research use only .