

产品名称: CD244 Mouse Monoclonal Antibody
产品货号: AMM82806



产品概述 (Summary)

产品名称 (Production Name)	CD244 Mouse Monoclonal Antibody
描述 (Description)	Mouse monoclonal Antibody
宿主 (Host)	Mouse
应用 (Application)	ICC,ELISA,FC
种属反应性 (Reactivity)	Human

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	Mouse IgG2b
克隆 (Clonality)	Monoclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	Purified antibody in PBS with 0.05% sodium azide
纯化方式 (Purification)	Affinity Purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	CD244
别名 (Alternative Names)	2B4; NAIL; Nmrk; NKR2B4; SLAMF4
基因 ID (Gene ID)	51744.0
蛋白 ID (SwissProt ID)	Q9BZW8.Purified recombinant fragment of human CD244 (AA: 22-229) expressed in E. Coli.

产品应用 (Application)

稀释比 (Dilution Ratio)	ICC 1:50-1:200,ELISA 1:5000-1:20000,FC 1:200-1:400
蛋白分子量 (Molecular Weight)	41.6kDa

研究背景 (Background)

This gene encodes a cell surface receptor expressed on natural killer (NK) cells (and some T cells) that mediate non-major


产品名称: CD244 Mouse Monoclonal Antibody
产品货号: AMM82806




histocompatibility complex (MHC) restricted killing. The interaction between NK-cell and target cells via this receptor is thought to modulate NK-cell cytolytic activity. Alternatively spliced transcript variants encoding different isoforms have been found for this gene.

研究领域 (Research Area)


图片 (Image Data)



Black line: Control Antigen (100 ng);Purple line: Antigen (10ng); Blue line: Antigen (50 ng); Red line:Antigen (100 ng)



Immunofluorescence analysis of HeLa cells using CD244 mouse mAb (green). Blue: DRAQ5 fluorescent DNA dye. Red: Actin filaments have been labeled with Alexa Fluor- 555 phalloidin. Secondary antibody from Fisher (Cat#: 35503)



产品名称: CD244 Mouse Monoclonal Antibody
产品货号: AMM82806



Flow cytometric analysis of THP-1 cells using CD244 mouse mAb (green) and negative control (red).

注意事项 (Note)

For research use only .