

产品名称: ITGB5 Mouse Monoclonal Antibody
产品货号: AMM82816



产品概述 (Summary)

产品名称 (Production Name)	ITGB5 Mouse Monoclonal Antibody
描述 (Description)	Mouse monoclonal Antibody
宿主 (Host)	Mouse
应用 (Application)	ICC,ELISA,FC
种属反应性 (Reactivity)	Human

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	Mouse IgG2a
克隆 (Clonality)	Monoclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	Purified antibody in PBS with 0.05% sodium azide
纯化方式 (Purification)	Affinity Purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	ITGB5
别名 (Alternative Names)	ITGB5
基因 ID (Gene ID)	3693.0
蛋白 ID (SwissProt ID)	P18084.Purified recombinant fragment of human ITGB5 (AA: 50-319) expressed in E. Coli.

产品应用 (Application)

稀释比 (Dilution Ratio)	ICC 1:50-1:200,ELISA 1:5000-1:20000,FC 1:200-1:400
蛋白分子量 (Molecular Weight)	88kDa

研究背景 (Background)

This gene encodes an integral outer mitochondrial membrane protein that blocks the apoptotic death of some cells such as

产品名称: ITGB5 Mouse Monoclonal Antibody
产品货号: AMM82816




lymphocytes. Constitutive expression of BCL2, such as in the case of translocation of BCL2 to Ig heavy chain locus, is thought to be the cause of follicular lymphoma. Alternative splicing results in multiple transcript variants.


研究领域 (Research Area)

PI3K-Akt signaling pathway


图片 (Image Data)



Black line: Control Antigen (100 ng);Purple line: Antigen (10ng); Blue line: Antigen (50 ng); Red line:Antigen (100 ng)



Immunofluorescence analysis of HeLa cells using ITGB5 mouse mAb (green). Blue: DRAQ5 fluorescent DNA dye. Red: Actin filaments have been labeled with Alexa Fluor- 555 phalloidin. Secondary antibody from Fisher (Cat#: 35503)



Flow cytometric analysis of Jurkat cells using ITGB5 mouse mAb (green) and negative control (red).

产品名称: ITGB5 Mouse Monoclonal Antibody
产品货号: AMM82816



注意事项 (Note)

For research use only .