

产品名称: ACADM Mouse Monoclonal Antibody
产品货号: AMM82819



产品概述 (Summary)

产品名称 (Production Name)	ACADM Mouse Monoclonal Antibody
描述 (Description)	Mouse monoclonal Antibody
宿主 (Host)	Mouse
应用 (Application)	WB,IHC,ICC,ELISA,FC
种属反应性 (Reactivity)	Human

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	Mouse IgG1
克隆 (Clonality)	Monoclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	Purified antibody in PBS with 0.05% sodium azide
纯化方式 (Purification)	Affinity Purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	ACADM
别名 (Alternative Names)	MCAD; ACAD1; MCADH
基因 ID (Gene ID)	34.0
蛋白 ID (SwissProt ID)	P11310.Purified recombinant fragment of human ACADM (AA: 26-185) expressed in E. Coli.

产品应用 (Application)

稀释比 (Dilution Ratio)	WB 1:500-1:2000,IHC 1:200-1:1000,ICC 1:200-1:1000,ELISA 1:5000-1:20000,FC 1:200-1:400
蛋白分子量 (Molecular Weight)	46.6KDa

研究背景 (Background)


产品名称: ACADM Mouse Monoclonal Antibody
产品货号: AMM82819




This gene encodes the medium-chain specific (C4 to C12 straight chain) acyl-Coenzyme A dehydrogenase. The homotetramer enzyme catalyzes the initial step of the mitochondrial fatty acid beta-oxidation pathway. Defects in this gene cause medium-chain acyl-CoA dehydrogenase deficiency, a disease characterized by hepatic dysfunction, fasting hypoglycemia, and encephalopathy, which can result in infantile death. Alternatively spliced transcript variants encoding different isoforms have been found for this gene.

研究领域 (Research Area)

图片 (Image Data)




Black line: Control Antigen (100 ng);Purple line: Antigen (10ng); Blue line: Antigen (50 ng); Red line:Antigen (100 ng)




Western blot analysis using ACADM mouse mAb against Hela (1), HepG2 (2), Jurkat (3), Raji (4) and K562 (5) cell lysate.


产品名称: ACADM Mouse Monoclonal Antibody
产品货号: AMM82819



Immunofluorescence analysis of Hela cells using ACADM mouse mAb (green). Blue: DRAQ5 fluorescent DNA dye. Red: Actin filaments have been labeled with Alexa Fluor- 555 phalloidin. Secondary antibody from Fisher (Cat#: 35503)



Flow cytometric analysis of HL-60 cells using ACADM mouse mAb (green) and negative control (red).



Immunohistochemical analysis of paraffin-embedded human cervical carcinoma tissues using ACADM mouse mAb with DAB staining.

注意事项 (Note)

For research use only .