

产品名称: MRGPRX2 Mouse Monoclonal Antibody
产品货号: AMM82842



产品概述 (Summary)

产品名称 (Production Name)	MRGPRX2 Mouse Monoclonal Antibody
描述 (Description)	Mouse monoclonal Antibody
宿主 (Host)	Mouse
应用 (Application)	ELISA,FC
种属反应性 (Reactivity)	Human

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	Mouse IgG1
克隆 (Clonality)	Monoclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	Purified antibody in PBS with 0.05% sodium azide
纯化方式 (Purification)	Affinity Purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	MRGPRX2
别名 (Alternative Names)	MGRG3; MRGX2
基因 ID (Gene ID)	117194.0
蛋白 ID (SwissProt ID)	Q96LB1.Purified recombinant fragment of human MRGPRX2 expressed in E. Coli.

产品应用 (Application)

稀释比 (Dilution Ratio)	ELISA 1:5000-1:20000,FC 1:200-1:400
蛋白分子量 (Molecular Weight)	37kDa

研究背景 (Background)

MRGPRX2 (MAS Related GPR Family Member X2) is a Protein Coding gene. Diseases associated with MRGPRX2 include


产品名称: MRGPRX2 Mouse Monoclonal Antibody
产品货号: AMM82842




Ceftazidime Allergy and Chronic Urticaria. Gene Ontology (GO) annotations related to this gene include G protein-coupled receptor activity and mast cell secretagogue receptor activity. An important paralog of this gene is MRGPRX1.

研究领域 (Research Area)

图片 (Image Data)



Black line: Control Antigen (100 ng);Purple line: Antigen (10ng); Blue line: Antigen (50 ng); Red line:Antigen (100 ng)



Flow cytometric analysis of MCF-7 cells using MRGPRX2 mouse mAb (green) and negative control (red).

注意事项 (Note)

For research use only .