

产品名称: DYRK4 Mouse Monoclonal Antibody
产品货号: AMM82862



产品概述 (Summary)

产品名称 (Production Name)	DYRK4 Mouse Monoclonal Antibody
描述 (Description)	Mouse monoclonal Antibody
宿主 (Host)	Mouse
应用 (Application)	IHC,ELISA
种属反应性 (Reactivity)	Human

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	Mouse IgG1
克隆 (Clonality)	Monoclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	Purified antibody in PBS with 0.05% sodium azide
纯化方式 (Purification)	Affinity Purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	DYRK4
别名 (Alternative Names)	DYRK4
基因 ID (Gene ID)	8798.0
蛋白 ID (SwissProt ID)	Q9NR20 .Purified recombinant fragment of human DYRK4 (AA: 303-515) expressed in E. Coli.

产品应用 (Application)

稀释比 (Dilution Ratio)	IHC 1:200-1:1000,ELISA 1:5000-1:20000
蛋白分子量 (Molecular Weight)	59.6kDa

研究背景 (Background)

This gene encodes an enzyme that belongs to a conserved family of serine/threonine protein kinases. Members of this dual


产品名称: DYRK4 Mouse Monoclonal Antibody
产品货号: AMM82862




specificity kinase family are thought to function in the regulation of cell differentiation and proliferation, survival, and in development. Alternate splicing results in multiple transcript variants. Additional alternatively spliced transcript variants of this gene have been described, but their full-length nature is not known.

研究领域 (Research Area)


图片 (Image Data)



Black line: Control Antigen (100 ng);Purple line: Antigen (10ng); Blue line: Antigen (50 ng); Red line:Antigen (100 ng)



Immunohistochemical analysis of paraffin-embedded human liver cancer tissues using DYRK4 mouse mAb with DAB staining.



产品名称: DYRK4 Mouse Monoclonal Antibody
产品货号: AMM82862



Immunohistochemical analysis of paraffin-embedded human pancreas tissues using DYRK4 mouse mAb with DAB staining.

注意事项 (Note)

For research use only .