

产品名称: TLT2 Mouse Monoclonal Antibody
产品货号: AMM82879



产品概述 (Summary)

产品名称 (Production Name)	TLT2 Mouse Monoclonal Antibody
描述 (Description)	Mouse monoclonal Antibody
宿主 (Host)	Mouse
应用 (Application)	IHC, ICC, ELISA, FC
种属反应性 (Reactivity)	Human

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	Mouse IgG2a
克隆 (Clonality)	Monoclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	Purified antibody in PBS with 0.05% sodium azide
纯化方式 (Purification)	Affinity Purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	TLT2
别名 (Alternative Names)	TLT2; TLT-2; C6orf76; dJ238O23.1
基因 ID (Gene ID)	79865.0
蛋白 ID (SwissProt ID)	Q5T2D2 .Purified recombinant fragment of human TLT2 (AA: extra 19-268) expressed in E. Coli.

产品应用 (Application)

稀释比 (Dilution Ratio)	IHC 1:200-1:1000, ICC 1:200-1:1000, ELISA 1:5000-1:20000, FC 1:200-1:400
蛋白分子量 (Molecular Weight)	35kDa

研究背景 (Background)

TREML2 is located in a gene cluster on chromosome 6 with the single Ig variable (IgV) domain activating receptors TREM1


产品名称: TLT2 Mouse Monoclonal Antibody
产品货号: AMM82879




(MIM 605085) and TREM2 (MIM 605086), but it has distinct structural and functional properties.

研究领域 (Research Area)


图片 (Image Data)



Black line: Control Antigen (100 ng);Purple line: Antigen (10ng); Blue line: Antigen (50 ng); Red line:Antigen (100 ng)




Immunofluorescence analysis of HeLa cells using TLT2 mouse mAb (green). Blue: DRAQ5 fluorescent DNA dye. Red: Actin filaments have been labeled with Alexa Fluor- 555 phalloidin. Secondary antibody from Fisher (Cat#: 35503)




Flow cytometric analysis of HEK293 cells using TLT2 mouse mAb (green) and negative control (red).

产品名称: TLT2 Mouse Monoclonal Antibody
产品货号: AMM82879



Immunohistochemical analysis of paraffin-embedded human cervical cancer tissues using TLT2 mouse mAb with DAB staining.



Immunohistochemical analysis of paraffin-embedded human rectum cancer tissues using TLT2 mouse mAb with DAB staining.

注意事项 (Note)

For research use only .