

产品名称: TYRO3 Mouse Monoclonal Antibody
产品货号: AMM82883



产品概述 (Summary)

产品名称 (Production Name)	TYRO3 Mouse Monoclonal Antibody
描述 (Description)	Mouse monoclonal Antibody
宿主 (Host)	Mouse
应用 (Application)	IHC,ELISA,FC
种属反应性 (Reactivity)	Human

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	Mouse IgG1
克隆 (Clonality)	Monoclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	Purified antibody in PBS with 0.05% sodium azide
纯化方式 (Purification)	Affinity Purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	TYRO3
别名 (Alternative Names)	BYK; Dtk; RSE; Rek; Sky; Tif; Etk-2
基因 ID (Gene ID)	7301.0
蛋白 ID (SwissProt ID)	Q06418.Purified recombinant fragment of human TYRO3 (AA: extra 230-429) expressed in mammalian.

产品应用 (Application)

稀释比 (Dilution Ratio)	IHC 1:200-1:1000,ELISA 1:5000-1:20000,FC 1:200-1:400
蛋白分子量 (Molecular Weight)	96.9kDa

研究背景 (Background)

The gene is part of a 3-member transmembrane receptor kinase receptor family with a processed pseudogene distal on


产品名称: TYRO3 Mouse Monoclonal Antibody
产品货号: AMM82883




chromosome 15. The encoded protein is activated by the products of the growth arrest-specific gene 6 and protein S genes and is involved in controlling cell survival and proliferation, spermatogenesis, immunoregulation and phagocytosis. The encoded protein has also been identified as a cell entry factor for Ebola and Marburg viruses.

研究领域 (Research Area)


图片 (Image Data)



Black line: Control Antigen (100 ng);Purple line: Antigen (10ng); Blue line: Antigen (50 ng); Red line:Antigen (100 ng)



Flow cytometric analysis of U937 cells using TYRO3 mouse mAb (green) and negative control (red).




Immunohistochemical analysis of paraffin-embedded human liver cancer tissues using TYRO3 mouse mAb with DAB

产品名称: TYRO3 Mouse Monoclonal Antibody
产品货号: AMM82883



staining.



Immunohistochemical analysis of paraffin-embedded human brain tissues using TYRO3 mouse mAb with DAB staining.

注意事项 (Note)

For research use only .