

产品名称: RBMXL3 Mouse Monoclonal Antibody
产品货号: AMM82891



产品概述 (Summary)

产品名称 (Production Name)	RBMXL3 Mouse Monoclonal Antibody
描述 (Description)	Mouse monoclonal Antibody
宿主 (Host)	Mouse
应用 (Application)	IHC, ICC, ELISA, FC
种属反应性 (Reactivity)	Human

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	Mouse IgG1
克隆 (Clonality)	Monoclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	Purified antibody in PBS with 0.05% sodium azide
纯化方式 (Purification)	Affinity Purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	RBMXL3
别名 (Alternative Names)	CXorf55
基因 ID (Gene ID)	139804.0
蛋白 ID (SwissProt ID)	Q8N7X1. Purified recombinant fragment of human RBMXL3(AA: 862-1067) expressed in E. Coli.

产品应用 (Application)

稀释比 (Dilution Ratio)	IHC 1:200-1:1000, ICC 1:50-1:200, ELISA 1:5000-1:20000, FC 1:200-1:400
蛋白分子量 (Molecular Weight)	114kDa

研究背景 (Background)

RBMXL3 (RBMX Like 3) is a Protein Coding gene. Diseases associated with RBMXL3 include Cardiomyopathy, Familial


产品名称: RBMXL3 Mouse Monoclonal Antibody
产品货号: AMM82891




Restrictive, 3. Among its related pathways are mRNA Splicing - Major Pathway. Gene Ontology (GO) annotations related to this gene include nucleic acid binding and nucleotide binding. An important paralog of this gene is RBMXL2.

研究领域 (Research Area)


图片 (Image Data)



Black line: Control Antigen (100 ng);Purple line: Antigen (10ng); Blue line: Antigen (50 ng); Red line:Antigen (100 ng)




Immunofluorescence analysis of HeLa cells using RBMXL3 mouse mAb (green). Blue: DRAQ5 fluorescent DNA dye. Red: Actin filaments have been labeled with Alexa Fluor- 555 phalloidin. Secondary antibody from Fisher (Cat#: 35503)




Flow cytometric analysis of Raji cells using RBMXL3 mouse mAb (green) and negative control (red).

产品名称: RBMXL3 Mouse Monoclonal Antibody
产品货号: AMM82891



Immunohistochemical analysis of paraffin-embedded human testis tissues using RBMXL3 mouse mAb with DAB staining.



Immunohistochemical analysis of paraffin-embedded human cerebellum tissues using RBMXL3 mouse mAb with DAB staining.

注意事项 (Note)

For research use only .