

产品名称: POU5F1 Mouse Monoclonal Antibody
产品货号: AMM82903



产品概述 (Summary)

产品名称 (Production Name)	POU5F1 Mouse Monoclonal Antibody
描述 (Description)	Mouse monoclonal Antibody
宿主 (Host)	Mouse
应用 (Application)	IHC,ELISA,FC
种属反应性 (Reactivity)	Human

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	Mouse IgG1
克隆 (Clonality)	Monoclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	Purified antibody in PBS with 0.05% sodium azide
纯化方式 (Purification)	Affinity Purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	POU5F1
别名 (Alternative Names)	OCT3; OCT4; OTF3; OTF4; OTF-3; Oct-3; Oct-4
基因 ID (Gene ID)	5460.0
蛋白 ID (SwissProt ID)	Q01860.Purified recombinant fragment of human (AA: 136-360) expressed in E. Coli.

产品应用 (Application)

稀释比 (Dilution Ratio)	IHC 1:200-1:1000,ELISA 1:5000-1:20000,FC 1:200-1:400
蛋白分子量 (Molecular Weight)	38.5kDa

研究背景 (Background)

This gene encodes a transcription factor containing a POU homeodomain that plays a key role in embryonic development


产品名称: POU5F1 Mouse Monoclonal Antibody
产品货号: AMM82903




and stem cell pluripotency. Aberrant expression of this gene in adult tissues is associated with tumorigenesis. This gene can participate in a translocation with the Ewing's sarcoma gene on chromosome 21, which also leads to tumor formation. Alternative splicing, as well as usage of alternative AUG and non-AUG translation initiation codons, results in multiple isoforms. One of the AUG start codons is polymorphic in human populations. Related pseudogenes have been identified on chromosomes 1, 3, 8, 10, and 12.

研究领域 (Research Area)

图片 (Image Data)




Black line: Control Antigen (100 ng);Purple line: Antigen (10ng); Blue line: Antigen (50 ng); Red line:Antigen (100 ng)




Flow cytometric analysis of HeLa cells using POU5F1 mouse mAb (green) and negative control (red).


产品名称: POU5F1 Mouse Monoclonal Antibody
产品货号: AMM82903



Flow cytometric analysis of K562 cells using POU5F1 mouse mAb (green) and negative control (red).



Immunohistochemical analysis of paraffin-embedded human colon cancer tissues using POU5F1 mouse mAb with DAB staining.



Immunohistochemical analysis of paraffin-embedded human rectum cancer tissues using POU5F1 mouse mAb with DAB staining.

注意事项 (Note)

For research use only .