

产品名称: NPM2 Mouse Monoclonal Antibody
产品货号: AMM82911



产品概述 (Summary)

产品名称 (Production Name)	NPM2 Mouse Monoclonal Antibody
描述 (Description)	Mouse monoclonal Antibody
宿主 (Host)	Mouse
应用 (Application)	IHC,ICC,ELISA
种属反应性 (Reactivity)	Human

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	Mouse IgG2b
克隆 (Clonality)	Monoclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	Purified antibody in PBS with 0.05% sodium azide
纯化方式 (Purification)	Affinity Purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	NPM2
别名 (Alternative Names)	NPM2
基因 ID (Gene ID)	10361.0
蛋白 ID (SwissProt ID)	Q86SE8 .Purified recombinant fragment of human NPM2 (AA: 1-214) expressed in E. Coli.

产品应用 (Application)

稀释比 (Dilution Ratio)	IHC 1:200-1:1000,ICC 1:100-1:400,ELISA 1:5000-1:20000
蛋白分子量 (Molecular Weight)	24.1kDa

研究背景 (Background)

NPM2 (Nucleophosmin/Nucleoplasmin 2) is a Protein Coding gene. Gene Ontology (GO) annotations related to this gene


产品名称: NPM2 Mouse Monoclonal Antibody
产品货号: AMM82911




include nucleic acid binding and enzyme binding. An important paralog of this gene is NPM1.

研究领域 (Research Area)


图片 (Image Data)



Black line: Control Antigen (100 ng);Purple line: Antigen (10ng); Blue line: Antigen (50 ng); Red line:Antigen (100 ng)




Immunofluorescence analysis of HeLa cells using NPM2 mouse mAb (green). Blue: DRAQ5 fluorescent DNA dye. Red: Actin filaments have been labeled with Alexa Fluor- 555 phalloidin. Secondary antibody from Fisher (Cat#: 35503)



Immunohistochemical analysis of paraffin-embedded human thyroid cancer tissues using NPM2 mouse mAb with DAB staining.

产品名称: NPM2 Mouse Monoclonal Antibody
产品货号: AMM82911



Immunohistochemical analysis of paraffin-embedded human brain tissues using NPM2 mouse mAb with DAB staining.

注意事项 (Note)

For research use only .