

产品名称: NPM2 Mouse Monoclonal Antibody
产品货号: AMM82912



产品概述 (Summary)

产品名称 (Production Name)	NPM2 Mouse Monoclonal Antibody
描述 (Description)	Mouse monoclonal Antibody
宿主 (Host)	Mouse
应用 (Application)	IHC,ELISA,FC
种属反应性 (Reactivity)	Human

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	Mouse IgG1
克隆 (Clonality)	Monoclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	Purified antibody in PBS with 0.05% sodium azide
纯化方式 (Purification)	Affinity Purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	NPM2
别名 (Alternative Names)	NPM2
基因 ID (Gene ID)	10361.0
蛋白 ID (SwissProt ID)	Q86SE8 .Purified recombinant fragment of human NPM2 (AA: 1-214) expressed in E. Coli.

产品应用 (Application)

稀释比 (Dilution Ratio)	IHC 1:200-1:1000,ELISA 1:5000-1:20000,FC 1:200-1:400
蛋白分子量 (Molecular Weight)	24.1kDa

研究背景 (Background)

NPM2 (Nucleophosmin/Nucleoplasmin 2) is a Protein Coding gene. Gene Ontology (GO) annotations related to this gene


产品名称: NPM2 Mouse Monoclonal Antibody
产品货号: AMM82912




include nucleic acid binding and enzyme binding. An important paralog of this gene is NPM1.

研究领域 (Research Area)


图片 (Image Data)



Black line: Control Antigen (100 ng);Purple line: Antigen (10ng); Blue line: Antigen (50 ng); Red line:Antigen (100 ng)



Flow cytometric analysis of HepG2 cells using NPM2 mouse mAb (green) and negative control (red).



Immunohistochemical analysis of paraffin-embedded human thyroid cancer tissues using NPM2 mouse mAb with DAB staining.

产品名称: NPM2 Mouse Monoclonal Antibody
产品货号: AMM82912



注意事项 (Note)

For research use only .