

产品名称: S100A4 Mouse Monoclonal Antibody
产品货号: AMM82951



产品概述 (Summary)

产品名称 (Production Name)	S100A4 Mouse Monoclonal Antibody
描述 (Description)	Mouse monoclonal Antibody
宿主 (Host)	Mouse
应用 (Application)	IHC,ELISA,FC
种属反应性 (Reactivity)	Human, Mouse, Rat

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	Mouse IgG1
克隆 (Clonality)	Monoclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	Purified antibody in PBS with 0.05% sodium azide
纯化方式 (Purification)	Affinity Purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	S100A4
别名 (Alternative Names)	42A; 18A2; CAPL; FSP1; MTS1; P9KA; PEL98
基因 ID (Gene ID)	6275.0
蛋白 ID (SwissProt ID)	P26447.Purified recombinant fragment of human S100A4 (AA: 2-101) expressed in E. Coli.

产品应用 (Application)

稀释比 (Dilution Ratio)	IHC 1:200-1:1000,ELISA 1:5000-1:20000,FC 1:200-1:400
蛋白分子量 (Molecular Weight)	11.7kDa

研究背景 (Background)

The protein encoded by this gene is a member of the S100 family of proteins containing 2 EF-hand calcium-binding motifs.


产品名称: S100A4 Mouse Monoclonal Antibody
产品货号: AMM82951




S100 proteins are localized in the cytoplasm and/or nucleus of a wide range of cells, and involved in the regulation of a number of cellular processes such as cell cycle progression and differentiation. S100 genes include at least 13 members which are located as a cluster on chromosome 1q21. This protein may function in motility, invasion, and tubulin polymerization. Chromosomal rearrangements and altered expression of this gene have been implicated in tumor metastasis. Multiple alternatively spliced variants, encoding the same protein, have been identified.

研究领域 (Research Area)

图片 (Image Data)




Black line: Control Antigen (100 ng);Purple line: Antigen (10ng); Blue line: Antigen (50 ng); Red line:Antigen (100 ng)




Flow cytometric analysis of Hela cells using S100A4 mouse mAb (green) and negative control (red).


产品名称: S100A4 Mouse Monoclonal Antibody
产品货号: AMM82951



Immunohistochemical analysis of paraffin-embedded human cervical cancer tissues using S100A4 mouse mAb with DAB staining.




Immunohistochemical analysis of paraffin-embedded human brain tissues using S100A4 mouse mAb with DAB staining.




Immunohistochemical analysis of paraffin-embedded mouse spleen tissues using S100A4 mouse mAb with DAB staining.


产品名称: S100A4 Mouse Monoclonal Antibody
产品货号: AMM82951




Immunohistochemical analysis of paraffin-embedded mouse kidney tissues using S100A4 mouse mAb with DAB staining.



Immunohistochemical analysis of paraffin-embedded rat spleen tissues using S100A4 mouse mAb with DAB staining.



Immunohistochemical analysis of paraffin-embedded rat kidney tissues using S100A4 mouse mAb with DAB staining.



产品名称: S100A4 Mouse Monoclonal Antibody
产品货号: AMM82951



Immunohistochemical analysis of paraffin-embedded rabbit brain tissues using S100A4 mouse mAb with DAB staining.

注意事项 (Note)

For research use only .