

产品名称: NOTCH4 Mouse Monoclonal Antibody
产品货号: AMM82959



产品概述 (Summary)

产品名称 (Production Name)	NOTCH4 Mouse Monoclonal Antibody
描述 (Description)	Mouse monoclonal Antibody
宿主 (Host)	Mouse
应用 (Application)	ICC,ELISA,FC
种属反应性 (Reactivity)	Human

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	Mouse IgG1
克隆 (Clonality)	Monoclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	Purified antibody in PBS with 0.05% sodium azide
纯化方式 (Purification)	Affinity Purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	NOTCH4
别名 (Alternative Names)	INT3
基因 ID (Gene ID)	4855.0
蛋白 ID (SwissProt ID)	Q99466.Purified recombinant fragment of human NOTCH4 (AA: extra 118-357) expressed in E. Coli.

产品应用 (Application)

稀释比 (Dilution Ratio)	ICC 1:50-1:200,ELISA 1:5000-1:20000,FC 1:200-1:400
蛋白分子量 (Molecular Weight)	209.6kDa

研究背景 (Background)

This gene encodes a member of the NOTCH family of proteins. Members of this Type I transmembrane protein family share

产品名称: NOTCH4 Mouse Monoclonal Antibody
产品货号: AMM82959




structural characteristics including an extracellular domain consisting of multiple epidermal growth factor-like (EGF) repeats, and an intracellular domain consisting of multiple different domain types. Notch signaling is an evolutionarily conserved intercellular signaling pathway that regulates interactions between physically adjacent cells through binding of Notch family receptors to their cognate ligands. The encoded preproprotein is proteolytically processed in the trans-Golgi network to generate two polypeptide chains that heterodimerize to form the mature cell-surface receptor. This receptor may play a role in vascular, renal and hepatic development. Mutations in this gene may be associated with schizophrenia. Alternative splicing results in multiple transcript variants, at least one of which encodes an isoform that is proteolytically processed.


研究领域 (Research Area)

Notch signaling pathway

图片 (Image Data)




Black line: Control Antigen (100 ng);Purple line: Antigen (10ng); Blue line: Antigen (50 ng); Red line:Antigen (100 ng)




Immunofluorescence analysis of HeLa cells using NOTCH4 mouse mAb (green). Blue: DRAQ5 fluorescent DNA dye. Red: Actin filaments have been labeled with Alexa Fluor- 555 phalloidin. Secondary antibody from Fisher (Cat#: 35503)

产品名称: NOTCH4 Mouse Monoclonal Antibody
产品货号: AMM82959



Flow cytometric analysis of Hela cells using NOTCH4 mouse mAb (green) and negative control (red).



Flow cytometric analysis of Jurkat cells using NOTCH4 mouse mAb (green) and negative control (red).

注意事项 (Note)

For research use only .