

产品名称: SMAD1 Mouse Monoclonal Antibody
产品货号: AMM82987



产品概述 (Summary)

产品名称 (Production Name)	SMAD1 Mouse Monoclonal Antibody
描述 (Description)	Mouse monoclonal Antibody
宿主 (Host)	Mouse
应用 (Application)	WB,ELISA
种属反应性 (Reactivity)	Human, Mouse, Rat

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	Mouse IgG1
克隆 (Clonality)	Monoclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	Purified antibody in PBS with 0.05% sodium azide
纯化方式 (Purification)	Affinity Purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	SMAD1
别名 (Alternative Names)	BSP1; JV41; BSP-1; JV4-1; MADH1; MADR1
基因 ID (Gene ID)	4086.0
蛋白 ID (SwissProt ID)	Q15797.Purified recombinant fragment of human SMAD1 (AA: 1-110) expressed in E. Coli.

产品应用 (Application)

稀释比 (Dilution Ratio)	WB 1:500-1:2000,ELISA 1:5000-1:20000
蛋白分子量 (Molecular Weight)	52.3kDa

研究背景 (Background)

The protein encoded by this gene belongs to the SMAD, a family of proteins similar to the gene products of the Drosophila

产品名称: SMAD1 Mouse Monoclonal Antibody
产品货号: AMM82987




gene 'mothers against decapentaplegic' (Mad) and the *C. elegans* gene Sma. SMAD proteins are signal transducers and transcriptional modulators that mediate multiple signaling pathways. This protein mediates the signals of the bone morphogenetic proteins (BMPs), which are involved in a range of biological activities including cell growth, apoptosis, morphogenesis, development and immune responses. In response to BMP ligands, this protein can be phosphorylated and activated by the BMP receptor kinase. The phosphorylated form of this protein forms a complex with SMAD4, which is important for its function in the transcription regulation. This protein is a target for SMAD-specific E3 ubiquitin ligases, such as SMURF1 and SMURF2, and undergoes ubiquitination and proteasome-mediated degradation. Alternatively spliced transcript variants encoding the same protein have been observed. [provided by RefSeq, Jul 2008]


研究领域 (Research Area)

TGF-beta signaling pathway

图片 (Image Data)



Black line: Control Antigen (100 ng);Purple line: Antigen (10ng); Blue line: Antigen (50 ng); Red line:Antigen (100 ng)



Western blot analysis using SMAD1 mouse mAb against Hela (1), HepG2 (2), MCF-7 (3), C2C12 (4), A431 (5), A549 (6), LNCap (7), NIH/3T3 (8) and PC-12 (9) cell lysate.

注意事项 (Note)

For research use only .

产品名称: SMAD1 Mouse Monoclonal Antibody
产品货号: AMM82987

