

产品名称: CD28 Mouse Monoclonal Antibody
产品货号: AMM82994



产品概述 (Summary)

产品名称 (Production Name)	CD28 Mouse Monoclonal Antibody
描述 (Description)	Mouse monoclonal Antibody
宿主 (Host)	Mouse
应用 (Application)	ICC,FC
种属反应性 (Reactivity)	Human,Mouse

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	Mouse IgG1
克隆 (Clonality)	Monoclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	Purified antibody in PBS with 0.05% sodium azide
纯化方式 (Purification)	Affinity Purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	CD28
别名 (Alternative Names)	Tp44
基因 ID (Gene ID)	940.0
蛋白 ID (SwissProt ID)	P10747.

产品应用 (Application)

稀释比 (Dilution Ratio)	ICC 1:25-1:50,FC 1:200-1:400
蛋白分子量 (Molecular Weight)	25kDa

研究背景 (Background)

The protein encoded by this gene is essential for T-cell proliferation and survival, cytokine production, and T-helper type-2 development. Several alternatively spliced transcript variants encoding different isoforms have been found for this gene.


产品名称: CD28 Mouse Monoclonal Antibody
产品货号: AMM82994



[provided by RefSeq, Jul 2011]

研究领域 (Research Area)

图片 (Image Data)



Confocal Immunofluorescence analysis of NIH3T3 cells using CD28 mouse mAb (green). Red: Actin filaments have been labeled with Alexa Fluor-555 phalloidin. Blue: DRAQ5 fluorescent DNA dye.

注意事项 (Note)

For research use only .