

产品名称: IL1A Mouse Monoclonal Antibody
产品货号: AMM83007



产品概述 (Summary)

产品名称 (Production Name)	IL1A Mouse Monoclonal Antibody
描述 (Description)	Mouse monoclonal Antibody
宿主 (Host)	Mouse
应用 (Application)	ELISA,FC
种属反应性 (Reactivity)	Human

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	Mouse IgG1
克隆 (Clonality)	Monoclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	Purified antibody in PBS with 0.05% sodium azide
纯化方式 (Purification)	Affinity Purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	IL1A
别名 (Alternative Names)	IL1; IL-1A; IL1F1; IL1-ALPHA; IL-1 alpha
基因 ID (Gene ID)	3552.0
蛋白 ID (SwissProt ID)	.Purified recombinant fragment of human IL1A (AA: extra(155-354)) expressed in E. Coli.

产品应用 (Application)

稀释比 (Dilution Ratio)	ELISA 1:5000-1:20000,FC 1:200-1:400
蛋白分子量 (Molecular Weight)	30.6kDa

研究背景 (Background)

The protein encoded by this gene is a member of the interleukin 1 cytokine family. This cytokine is a pleiotropic cytokine

产品名称: IL1A Mouse Monoclonal Antibody
产品货号: AMM83007




involved in various immune responses, inflammatory processes, and hematopoiesis. This cytokine is produced by monocytes and macrophages as a proprotein, which is proteolytically processed and released in response to cell injury, and thus induces apoptosis. This gene and eight other interleukin 1 family genes form a cytokine gene cluster on chromosome 2. It has been suggested that the polymorphism of these genes is associated with rheumatoid arthritis and Alzheimer's disease.


研究领域 (Research Area)

TGF-beta signaling pathway, MAPK signaling pathway

图片 (Image Data)




Black line: Control Antigen (100 ng); Purple line: Antigen (10ng); Blue line: Antigen (50 ng); Red line: Antigen (100 ng)




Flow cytometric analysis of Jurkat cells using IL1A mouse mAb (green) and negative control (red).

产品名称: IL1A Mouse Monoclonal Antibody
产品货号: AMM83007



Flow cytometric analysis of MOLT4 cells using IL1A mouse mAb (green) and negative control (red).



Flow cytometric analysis of THP-1 cells using IL1A mouse mAb (green) and negative control (red).

注意事项 (Note)

For research use only .