

产品名称: GNL3 Mouse Monoclonal Antibody
产品货号: AMM83041



产品概述 (Summary)

产品名称 (Production Name)	GNL3 Mouse Monoclonal Antibody
描述 (Description)	Mouse monoclonal Antibody
宿主 (Host)	Mouse
应用 (Application)	WB,ELISA,FC
种属反应性 (Reactivity)	Human,Mouse

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	Mouse IgG1
克隆 (Clonality)	Monoclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	Purified antibody in PBS with 0.05% sodium azide
纯化方式 (Purification)	Affinity Purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	GNL3
别名 (Alternative Names)	NS; E2IG3; NNP47; C77032
基因 ID (Gene ID)	26354.0
蛋白 ID (SwissProt ID)	Q9BVP2.Purified recombinant fragment of human GNL3 (AA: 1-226) expressed in E. Coli.

产品应用 (Application)

稀释比 (Dilution Ratio)	WB 1:500-1:2000,ELISA 1:5000-1:20000,FC 1:200-1:400
蛋白分子量 (Molecular Weight)	62kDa

研究背景 (Background)

The protein encoded by this gene may interact with p53 and may be involved in tumorigenesis. The encoded protein also


产品名称: GNL3 Mouse Monoclonal Antibody
产品货号: AMM83041




appears to be important for stem cell proliferation. This protein is found in both the nucleus and nucleolus. Three transcript variants encoding two different isoforms have been found for this gene.

研究领域 (Research Area)


图片 (Image Data)



Black line: Control Antigen (100 ng); Purple line: Antigen(10ng); Blue line: Antigen (50 ng); Red line: Antigen (100 ng);



Western blot analysis using GNL3 mouse mAb against NIH3T3 (1) and PC-3 (2) cell lysate.



Flow cytometric analysis of Jurkat cells using GNL3 mouse mAb (green) and negative control (purple).

产品名称: GNL3 Mouse Monoclonal Antibody
产品货号: AMM83041



注意事项 (Note)

For research use only .