

产品名称: TRAFD1 Mouse Monoclonal Antibody
产品货号: AMM83043



产品概述 (Summary)

产品名称 (Production Name)	TRAFD1 Mouse Monoclonal Antibody
描述 (Description)	Mouse monoclonal Antibody
宿主 (Host)	Mouse
应用 (Application)	ICC,ELISA,FC
种属反应性 (Reactivity)	Human

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	Mouse IgG1
克隆 (Clonality)	Monoclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	Purified antibody in PBS with 0.05% sodium azide
纯化方式 (Purification)	Affinity Purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	TRAFD1
别名 (Alternative Names)	FLN29
基因 ID (Gene ID)	10906.0
蛋白 ID (SwissProt ID)	O14545 .Purified recombinant fragment of human TRAFD1 (AA: 401-582) expressed in E. Coli.

产品应用 (Application)

稀释比 (Dilution Ratio)	ICC 1:100-1:200,ELISA 1:5000-1:20000,FC 1:200-1:400
蛋白分子量 (Molecular Weight)	64.8kDa

研究背景 (Background)

The innate immune system confers host defense against viral and microbial infection, and TRAFD1 is a negative feedback


产品名称: TRAFD1 Mouse Monoclonal Antibody
产品货号: AMM83043




regulator that controls excessive immune responses

研究领域 (Research Area)


图片 (Image Data)



Black line: Control Antigen (100 ng); Purple line: Antigen(10ng); Blue line: Antigen (50 ng); Red line: Antigen (100 ng);




Western blot analysis using TRAFD1 mouse mAb against HEK293 (1), Raji (2), and Jurkat (3) cell lysate.



Immunofluorescence analysis of HepG2 cells using TRAFD1 mouse mAb (green). Blue: DRAQ5 fluorescent DNA dye.

产品名称: TRAFD1 Mouse Monoclonal Antibody
产品货号: AMM83043



Flow cytometric analysis of HeLa cells using TRAFD1 mouse mAb (green) and negative control (purple).

注意事项 (Note)

For research use only .