

产品名称: PDE1B Mouse Monoclonal Antibody
产品货号: AMM83049



产品概述 (Summary)

产品名称 (Production Name)	PDE1B Mouse Monoclonal Antibody
描述 (Description)	Mouse monoclonal Antibody
宿主 (Host)	Mouse
应用 (Application)	ELISA,FC
种属反应性 (Reactivity)	Human,Mouse,Monkey

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	Mouse IgG1
克隆 (Clonality)	Monoclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	Purified antibody in PBS with 0.05% sodium azide
纯化方式 (Purification)	Affinity Purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	PDE1B
别名 (Alternative Names)	PDE1B1; PDES1B
基因 ID (Gene ID)	5153.0
蛋白 ID (SwissProt ID)	Q01064.Purified recombinant fragment of human PDE1B (AA: 370-536) expressed in E. Coli.

产品应用 (Application)

稀释比 (Dilution Ratio)	ELISA 1:5000-1:20000,FC 1:200-1:400
蛋白分子量 (Molecular Weight)	61.4kDa

研究背景 (Background)

The protein encoded by this gene belongs to the cyclic nucleotide phosphodiesterase (PDE) family, and PDE1 subfamily.


产品名称: PDE1B Mouse Monoclonal Antibody
产品货号: AMM83049




Members of the PDE1 family are calmodulin-dependent PDEs that are stimulated by a calcium-calmodulin complex. This PDE has dual-specificity for the second messengers, cAMP and cGMP, with a preference for cGMP as a substrate. cAMP and cGMP function as key regulators of many important physiological processes. Alternatively spliced transcript variants encoding different isoforms have been described for this gene.

研究领域 (Research Area)

图片 (Image Data)




Black line: Control Antigen (100 ng); Purple line: Antigen(10ng); Blue line: Antigen (50 ng); Red line: Antigen (100 ng);



Western blot analysis using PDE1B mouse mAb against A549 (1), SK-MES-1 (2), PC-12 (4), 3T3L1 (5) cell lysate and Mouse brain (3) tissue lysate.

产品名称: PDE1B Mouse Monoclonal Antibody
产品货号: AMM83049



Flow cytometric analysis of A549 cells using PDE1B mouse mAb (green) and negative control (red).

注意事项 (Note)

For research use only .