

产品名称: PIWIL4 Mouse Monoclonal Antibody
产品货号: AMM83056



产品概述 (Summary)

产品名称 (Production Name)	PIWIL4 Mouse Monoclonal Antibody
描述 (Description)	Mouse monoclonal Antibody
宿主 (Host)	Mouse
应用 (Application)	WB,IHC,ELISA,FC
种属反应性 (Reactivity)	Human

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	Mouse IgG1
克隆 (Clonality)	Monoclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	PBS containing 0.03% sodium azide.
纯化方式 (Purification)	Affinity Purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	PIWIL4
别名 (Alternative Names)	HIWI2; MIWI2
基因 ID (Gene ID)	143689.0
蛋白 ID (SwissProt ID)	Q7Z3Z4.Purified recombinant fragment of human PIWIL4 (AA: 304-434) expressed in E. Coli.

产品应用 (Application)

稀释比 (Dilution Ratio)	WB 1:500-1:2000,IHC 1:200-1:1000,ELISA 1:5000-1:20000,FC 1:200-1:400
蛋白分子量 (Molecular Weight)	96.6kDa

研究背景 (Background)

PIWIL4 belongs to the Argonaute family of proteins, which function in development and maintenance of germline stem


产品名称: PIWIL4 Mouse Monoclonal Antibody
产品货号: AMM83056



cells

研究领域 (Research Area)

图片 (Image Data)



Immunohistochemical analysis of paraffin-embedded human rectum cancer tissues using PIWIL4 mouse mAb with DAB staining.

注意事项 (Note)

For research use only .