

产品名称: DNAL4 Mouse Monoclonal Antibody
产品货号: AMM83071



产品概述 (Summary)

产品名称 (Production Name)	DNAL4 Mouse Monoclonal Antibody
描述 (Description)	Mouse monoclonal Antibody
宿主 (Host)	Mouse
应用 (Application)	IHC,ELISA
种属反应性 (Reactivity)	Human

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	Mouse IgG1
克隆 (Clonality)	Monoclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	Purified antibody in PBS with 0.05% sodium azide
纯化方式 (Purification)	Affinity Purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	DNAL4
别名 (Alternative Names)	PIG27
基因 ID (Gene ID)	10126.0
蛋白 ID (SwissProt ID)	O96015.Purified recombinant fragment of human DNAL4 (AA: 1-105) expressed in E. Coli.

产品应用 (Application)

稀释比 (Dilution Ratio)	IHC 1:200-1:1000,ELISA 1:5000-1:20000
蛋白分子量 (Molecular Weight)	12kDa

研究背景 (Background)


DNAL4 is a component of the dynein motor complex

产品名称: DNAL4 Mouse Monoclonal Antibody
产品货号: AMM83071



研究领域 (Research Area)

图片 (Image Data)



Red: Control Antigen (100ng);

注意事项 (Note)

For research use only .