

产品名称: PAX3 Mouse Monoclonal Antibody
产品货号: AMM83076



产品概述 (Summary)

产品名称 (Production Name)	PAX3 Mouse Monoclonal Antibody
描述 (Description)	Mouse monoclonal Antibody
宿主 (Host)	Mouse
应用 (Application)	WB,ELISA
种属反应性 (Reactivity)	Human,Mouse

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	Mouse IgG1
克隆 (Clonality)	Monoclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	PBS containing 0.03% sodium azide.
纯化方式 (Purification)	Affinity Purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	PAX3
别名 (Alternative Names)	WS1; WS3; CDHS; HUP2
基因 ID (Gene ID)	5077.0
蛋白 ID (SwissProt ID)	P23760.Purified recombinant fragment of human PAX3 (AA: 142-203) expressed in E. Coli.

产品应用 (Application)

稀释比 (Dilution Ratio)	WB 1:500-1:2000,ELISA 1:5000-1:20000
蛋白分子量 (Molecular Weight)	53kDa

研究背景 (Background)

This gene is a member of the paired box (PAX) family of transcription factors. Members of the PAX family typically contain


产品名称: PAX3 Mouse Monoclonal Antibody
产品货号: AMM83076




a paired box domain and a paired-type homeodomain. These genes play critical roles during fetal development. Mutations in paired box gene 3 are associated with Waardenburg syndrome, craniofacial-deafness-hand syndrome, and alveolar rhabdomyosarcoma. The translocation t(2;13)(q35;q14), which represents a fusion between PAX3 and the forkhead gene, is a frequent finding in alveolar rhabdomyosarcoma. Alternative splicing results in transcripts encoding isoforms with different C-termini.

研究领域 (Research Area)

图片 (Image Data)



Black line: Control Antigen (100 ng); Purple line: Antigen(10ng); Blue line: Antigen (50 ng); Red line: Antigen (100 ng);



Western blot analysis using PAX3 mouse mAb against A549 (1), A431 (2), Jurkat (3), Rat spleen (4) and Mouse liver (5) cell lysate.

注意事项 (Note)

For research use only .