

产品名称: TP53BP1 Mouse Monoclonal Antibody
产品货号: AMM83080



产品概述 (Summary)

产品名称 (Production Name)	TP53BP1 Mouse Monoclonal Antibody
描述 (Description)	Mouse monoclonal Antibody
宿主 (Host)	Mouse
应用 (Application)	WB,IHC,ELISA,FC
种属反应性 (Reactivity)	Human

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	Mouse IgG1
克隆 (Clonality)	Monoclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	PBS containing 0.03% sodium azide.
纯化方式 (Purification)	Affinity Purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	TP53BP1
别名 (Alternative Names)	p202; 53BP1
基因 ID (Gene ID)	7158.0
蛋白 ID (SwissProt ID)	Q12888.Purified recombinant fragment of human TP53BP1 (AA: 574-773) expressed in E. Coli.

产品应用 (Application)

稀释比 (Dilution Ratio)	WB 1:500-1:2000,IHC 1:200-1:1000,ELISA 1:5000-1:20000,FC 1:200-1:400
蛋白分子量 (Molecular Weight)	213.6kDa

研究背景 (Background)

T protein p53 binding protein 1 may have a role in checkpoint signaling during mitosis,enhance TP53-mediated


产品名称: TP53BP1 Mouse Monoclonal Antibody
产品货号: AMM83080




transcriptional activation and play a role in the response to DNA damage.

研究领域 (Research Area)


图片 (Image Data)



Black line: Control Antigen (100 ng); Purple line: Antigen(10ng); Blue line: Antigen (50 ng); Red line: Antigen (100 ng);




Flow cytometric analysis of HepG2 cells using TP53BP1 mouse mAb (green) and negative control (purple).



Immunohistochemical analysis of paraffin-embedded human colon cancer tissues using TP53BP1 mouse mAb with DAB staining.

产品名称: TP53BP1 Mouse Monoclonal Antibody
产品货号: AMM83080



Immunohistochemical analysis of paraffin-embedded human endometrial cancer tissues using TP53BP1 mouse mAb with DAB staining.

注意事项 (Note)

For research use only .