

产品名称: Human IgA Mouse Monoclonal Antibody
产品货号: AMM83086



产品概述 (Summary)

产品名称 (Production Name)	Human IgA Mouse Monoclonal Antibody
描述 (Description)	Mouse monoclonal Antibody
宿主 (Host)	Mouse
应用 (Application)	WB,IHC
种属反应性 (Reactivity)	Human

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	Mouse IgG1
克隆 (Clonality)	Monoclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	Purified antibody in PBS with 0.05% sodium azide
纯化方式 (Purification)	Affinity Purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	Human IgA
别名 (Alternative Names)	Human IgA
基因 ID (Gene ID)	
蛋白 ID (SwissProt ID)	P01877.Recombinant human IgA2

产品应用 (Application)

稀释比 (Dilution Ratio)	WB 1:500-1:1000,IHC 1:500-1:1000
蛋白分子量 (Molecular Weight)	36.6kDa


研究背景 (Background)

产品名称: Human IgA Mouse Monoclonal Antibody
产品货号: AMM83086



研究领域 (Research Area)

图片 (Image Data)



Microtiter wells were coated with human immunoglobulins with various isotypes. AMM83086 was used as the primary antibody and peroxidase conjugated goat anti-mouse IgG was used as the secondary antibody.

注意事项 (Note)

For research use only .