

产品名称: GBP6 Mouse Monoclonal Antibody
产品货号: AMM83092



产品概述 (Summary)

产品名称 (Production Name)	GBP6 Mouse Monoclonal Antibody
描述 (Description)	Mouse monoclonal Antibody
宿主 (Host)	Mouse
应用 (Application)	WB,IHC,IP
种属反应性 (Reactivity)	Human

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	Mouse IgG2a
克隆 (Clonality)	Monoclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	Purified antibody in PBS with 0.05% sodium azide
纯化方式 (Purification)	Affinity Purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	GBP6
别名 (Alternative Names)	GBP6
基因 ID (Gene ID)	
蛋白 ID (SwissProt ID)	Q6ZN66.Purified recombinant fragment of human GBP6

产品应用 (Application)

稀释比 (Dilution Ratio)	WB 1:500-1:1000,IHC 1:50-1:200,IP 1:100-1:200
蛋白分子量 (Molecular Weight)	72kDa


研究背景 (Background)

产品名称: GBP6 Mouse Monoclonal Antibody
产品货号: AMM83092




研究领域 (Research Area)


图片 (Image Data)



Western blot analysis using GBP6 mouse mAb against tonsil tissue lysate



Lane 1: Tonsil tissue lysate ; Lane 2: GBP6 immunoprecipitated from tonsil tissue lysate by AMM83092; Lane 3: The same as Lane 2 but use IgG isotype control antibody ;Result: AMM83092 can immunoprecipitate GBP6;



IHC-P analysis of esophagus tissue by GBP6 antibody (AMM83092)

注意事项 (Note)

产品名称: GBP6 Mouse Monoclonal Antibody
产品货号: AMM83092



For research use only .