

产品名称: SNAP25 Mouse Monoclonal Antibody
产品货号: AMM83095



产品概述 (Summary)

产品名称 (Production Name)	SNAP25 Mouse Monoclonal Antibody
描述 (Description)	Mouse monoclonal Antibody
宿主 (Host)	Mouse
应用 (Application)	WB,IHC,IP
种属反应性 (Reactivity)	Human

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	Mouse IgG1
克隆 (Clonality)	Monoclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	Purified antibody in PBS with 0.05% sodium azide
纯化方式 (Purification)	Affinity Purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	SNAP25
别名 (Alternative Names)	SNAP25
基因 ID (Gene ID)	
蛋白 ID (SwissProt ID)	P60880.Purified recombinant fragment of human SNAP25

产品应用 (Application)

稀释比 (Dilution Ratio)	WB 1:500-1:1000,IHC 1:50-1:250,IP 1:100-1:200
蛋白分子量 (Molecular Weight)	Calculated MW: 23 kDa; Observed MW: 23 kDa


研究背景 (Background)

产品名称: SNAP25 Mouse Monoclonal Antibody
产品货号: AMM83095




研究领域 (Research Area)


图片 (Image Data)



Lane 1: 15 µg of SH-SY5Y lysate; Lane 2: 15 µg of rat brain tissue lysate; Result: AMM83095 can detect SNAP25 by Western



Lane 1: Rat brain tissue lysate; Lane 2: SNAP25 immunoprecipitated from rat brain tissue lysate by AMM83095; Lane 3: The same as Lane 2 but use IgG isotype control antibody; Result: AMM83095 can immunoprecipitate SNAP25;



IHC-P analysis of cerebral cortex tissue by anti-SNAP25 antibody (AMM83095)

产品名称: SNAP25 Mouse Monoclonal Antibody
产品货号: AMM83095



注意事项 (Note)

For research use only .