

产品名称: DNAJB6 Mouse Monoclonal Antibody
产品货号: AMM83099



产品概述 (Summary)

产品名称 (Production Name)	DNAJB6 Mouse Monoclonal Antibody
描述 (Description)	Mouse monoclonal Antibody
宿主 (Host)	Mouse
应用 (Application)	WB,IHC,IP
种属反应性 (Reactivity)	Human

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	Mouse IgG2a
克隆 (Clonality)	Monoclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	Purified antibody in PBS with 0.05% sodium azide
纯化方式 (Purification)	Affinity Purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	DNAJB6
别名 (Alternative Names)	DNAJB6
基因 ID (Gene ID)	
蛋白 ID (SwissProt ID)	O75190.Purified recombinant fragment of human DNAJB6

产品应用 (Application)

稀释比 (Dilution Ratio)	WB 1:500-1:1000,IHC 1:100-1:200,IP 1:100-1:200
蛋白分子量 (Molecular Weight)	Calculated MW: 36 kDa; Observed MW: 36 kDa


研究背景 (Background)

产品名称: DNAJB6 Mouse Monoclonal Antibody
产品货号: AMM83099




研究领域 (Research Area)


图片 (Image Data)



Western blot analysis using DNAJB6 mouse mAb against SH-SY5Y lysate



Lane 1: SH-SY5Y lysate ; Lane 2: DNAJB6 immunoprecipitated from SH-SY5Y lysate by AMM83099 ; Lane 3: The same as Lane 2 but use IgG isotype control antibody ;Result: AMM83099 can immunoprecipitate DNAJB6;



IHC-P analysis of testis tissue by DNAJB6 antibody (AMM83099)

注意事项 (Note)

产品名称: DNAJB6 Mouse Monoclonal Antibody
产品货号: AMM83099



For research use only .