

产品名称: EXOSC10 Mouse Monoclonal Antibody
产品货号: AMM83106



产品概述 (Summary)

产品名称 (Production Name)	EXOSC10 Mouse Monoclonal Antibody
描述 (Description)	Mouse monoclonal Antibody
宿主 (Host)	Mouse
应用 (Application)	WB,IP
种属反应性 (Reactivity)	Human

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	Mouse IgG2a
克隆 (Clonality)	Monoclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	Purified antibody in PBS with 0.05% sodium azide
纯化方式 (Purification)	Affinity Purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	EXOSC10
别名 (Alternative Names)	EXOSC10
基因 ID (Gene ID)	
蛋白 ID (SwissProt ID)	Q01780.Purified recombinant fragment of human EXOSC10

产品应用 (Application)

稀释比 (Dilution Ratio)	WB 1:500-1:1000,IP 1:100-1:200
蛋白分子量 (Molecular Weight)	100.8kDa


研究背景 (Background)

产品名称: EXOSC10 Mouse Monoclonal Antibody
产品货号: AMM83106




研究领域 (Research Area)

图片 (Image Data)



Western blot analysis using EXOSC10 mouse mAb against HEK-293 lysate



Lane 1: HEK-293 lysate ; Lane 2: EXOSC10 immunoprecipitated from HEK-293 lysate by AMM83106 ; Lane 3: The same as Lane 2 but use IgG isotype control antibody ;Result: AMM83106 can immunoprecipitate EXOSC10;

注意事项 (Note)

For research use only .