

产品名称: ENO1 Mouse Monoclonal Antibody
产品货号: AMM83109



产品概述 (Summary)

产品名称 (Production Name)	ENO1 Mouse Monoclonal Antibody
描述 (Description)	Mouse monoclonal Antibody
宿主 (Host)	Mouse
应用 (Application)	WB,IHC,IP
种属反应性 (Reactivity)	Human

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	Mouse IgG1
克隆 (Clonality)	Monoclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	Purified antibody in PBS with 0.05% sodium azide
纯化方式 (Purification)	Affinity Purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	ENO1
别名 (Alternative Names)	ENO1
基因 ID (Gene ID)	
蛋白 ID (SwissProt ID)	P06733.Purified recombinant fragment of human ENO1

产品应用 (Application)

稀释比 (Dilution Ratio)	WB 1:500-1:1000,IHC 1:200-1:500,IP 1:100-1:200
蛋白分子量 (Molecular Weight)	Calculated MW: 47 kDa; Observed MW: 47 kDa


研究背景 (Background)

产品名称: ENO1 Mouse Monoclonal Antibody
产品货号: AMM83109




研究领域 (Research Area)


图片 (Image Data)



Lane 1: 15 μ g of SH-SY5Y lysate; Lane 2: 15 μ g of U-251MG lysate; Lane 3: 15 μ g of MCF-7 lysate; Lane 4: 15 μ g of A-431 lysate; Lane 5: 15 μ g of HeLa lysate; Lane 6: 15 μ g of THP-1 lysate; Lane 7: 15 μ g of Jurkat lysate; Lane 8: 15 μ g of Ramos lysate; Lane 9: 15 μ g of human plasma; Result: AMM83109 can detect human ENO1 by Western blotting.



Lane 1: 15 μ g of MCF-7 lysate; Lane 2: ENO1 immunoprecipitated from MCF-7 lysate by AMM83109; Lane 3: The same as Lane 2 but use IgG isotype control antibody; Result: AMM83109 can immunoprecipitate ENO1;



IHC-P analysis of colon tissue by ENO1 antibody (AMM83109)

产品名称: ENO1 Mouse Monoclonal Antibody
产品货号: AMM83109



注意事项 (Note)

For research use only .