

产品名称: CARD9 Mouse Monoclonal Antibody
产品货号: AMM83115



产品概述 (Summary)

产品名称 (Production Name)	CARD9 Mouse Monoclonal Antibody
描述 (Description)	Mouse monoclonal Antibody
宿主 (Host)	Mouse
应用 (Application)	WB,IP
种属反应性 (Reactivity)	Human

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	Mouse IgG1
克隆 (Clonality)	Monoclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	Purified antibody in PBS with 0.05% sodium azide
纯化方式 (Purification)	Affinity Purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	CARD9
别名 (Alternative Names)	CARD9
基因 ID (Gene ID)	
蛋白 ID (SwissProt ID)	Q9H257.Purified recombinant fragment of human CARD9

产品应用 (Application)

稀释比 (Dilution Ratio)	WB 1:500-1:1000,IP 1:100-1:200
蛋白分子量 (Molecular Weight)	Calculated MW: 62 kDa; Observed MW: 62 kDa


研究背景 (Background)

产品名称: CARD9 Mouse Monoclonal Antibody
产品货号: AMM83115




研究领域 (Research Area)

图片 (Image Data)



Western blot analysis using CARD9 mouse mAb against THP-1 lysate



Lane 1: THP-1 lysate ; Lane 2: CARD9 immunoprecipitated from THP-1 lysate by AMM83115 ; Lane 3: The same as Lane 2 but use IgG isotype control antibody ;Result: AMM83115 can immunoprecipitate CARD9;

注意事项 (Note)

For research use only .