

产品名称: LGALS2 Mouse Monoclonal Antibody
产品货号: AMM83126



产品概述 (Summary)

产品名称 (Production Name)	LGALS2 Mouse Monoclonal Antibody
描述 (Description)	Mouse monoclonal Antibody
宿主 (Host)	Mouse
应用 (Application)	WB,IHC,IP
种属反应性 (Reactivity)	Human

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	Mouse IgG1
克隆 (Clonality)	Monoclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	Purified antibody in PBS with 0.05% sodium azide
纯化方式 (Purification)	Affinity Purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	LGALS2
别名 (Alternative Names)	LGALS2
基因 ID (Gene ID)	
蛋白 ID (SwissProt ID)	P05162.Purified recombinant fragment of human LGALS2

产品应用 (Application)

稀释比 (Dilution Ratio)	WB 1:500-1:1000,IHC 1:100-1:200,IP 1:100-1:200
蛋白分子量 (Molecular Weight)	14.6kDa


研究背景 (Background)

产品名称: LGALS2 Mouse Monoclonal Antibody
产品货号: AMM83126




研究领域 (Research Area)


图片 (Image Data)



Western blot analysis using LGALS2 mouse mAb against kidney tissue lysate



Lane 1: Kidney tissue lysate ; Lane 2: LGALS2 immunoprecipitated from kidney tissue lysate by AMM83126 ; Lane 3: The same as Lane 2 but use IgG isotype control antibody ;Result: AMM83126 can immunoprecipitate LGALS2;



IHC-P analysis of gallbladder tissue by LGALS2 antibody (AMM83126)

注意事项 (Note)

产品名称: LGALS2 Mouse Monoclonal Antibody
产品货号: AMM83126



For research use only .