

产品名称: DGKZ Mouse Monoclonal Antibody
产品货号: AMM83130



产品概述 (Summary)

产品名称 (Production Name)	DGKZ Mouse Monoclonal Antibody
描述 (Description)	Mouse monoclonal Antibody
宿主 (Host)	Mouse
应用 (Application)	WB,IHC,IP
种属反应性 (Reactivity)	Human, Mouse

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	Mouse IgG2b
克隆 (Clonality)	Monoclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	Purified antibody in PBS with 0.05% sodium azide
纯化方式 (Purification)	Affinity Purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	DGKZ
别名 (Alternative Names)	DGKZ
基因 ID (Gene ID)	
蛋白 ID (SwissProt ID)	Q13574-1.Purified recombinant fragment of human DGKZ

产品应用 (Application)

稀释比 (Dilution Ratio)	WB 1:500-1:1000,IHC 1:100-1:200,IP 1:100-1:200
蛋白分子量 (Molecular Weight)	130kDa


研究背景 (Background)

产品名称: DGKZ Mouse Monoclonal Antibody
产品货号: AMM83130




研究领域 (Research Area)


图片 (Image Data)



Western blot analysis using DGKZ mouse mAb against mouse brain lysate



Lane 1: Mouse brain tissue lysate ; Lane 2: DGKZ immunoprecipitated from mouse brain tissue lysate by AMM83130 ; Lane 3: The same as Lane 2 but use IgG isotype control antibody ;Result: AMM83130 can immunoprecipitate DGKZ;



IHC-P analysis of cerebellum tissue by DGKZ antibody (AMM83130)

注意事项 (Note)

产品名称: DGKZ Mouse Monoclonal Antibody
产品货号: AMM83130



For research use only .