

产品名称: HSPA1A Mouse Monoclonal Antibody
产品货号: AMM83131



产品概述 (Summary)

产品名称 (Production Name)	HSPA1A Mouse Monoclonal Antibody
描述 (Description)	Mouse monoclonal Antibody
宿主 (Host)	Mouse
应用 (Application)	WB,IHC,IP
种属反应性 (Reactivity)	Human

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	Mouse IgG1
克隆 (Clonality)	Monoclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	Purified antibody in PBS with 0.05% sodium azide
纯化方式 (Purification)	Affinity Purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	HSPA1A
别名 (Alternative Names)	HSPA1A
基因 ID (Gene ID)	
蛋白 ID (SwissProt ID)	P0DMV8.Purified recombinant fragment of human HSPA1A

产品应用 (Application)

稀释比 (Dilution Ratio)	WB 1:500-1:1000,IHC 1:100-1:200,IP 1:100-1:200
蛋白分子量 (Molecular Weight)	Calculated MW: 70 kDa; Observed MW: 70 kDa


研究背景 (Background)

产品名称: HSPA1A Mouse Monoclonal Antibody
产品货号: AMM83131




研究领域 (Research Area)


图片 (Image Data)



Western blot analysis using HSPA1A mouse mAb against A549 lysate



Lane 1: A549 lysate ; Lane 2: HSPA1A immunoprecipitated from A549 lysate by AMM83131 ; Lane 3: The same as Lane 2 but use IgG isotype control antibody ;Result: AMM83131 can immunoprecipitate HSPA1A;



IHC-P analysis of testis tissue by HSPA1A antibody (AMM83131)

注意事项 (Note)

产品名称: HSPA1A Mouse Monoclonal Antibody
产品货号: AMM83131



For research use only .