

产品名称: C7orf50 Mouse Monoclonal Antibody
产品货号: AMM83133



产品概述 (Summary)

产品名称 (Production Name)	C7orf50 Mouse Monoclonal Antibody
描述 (Description)	Mouse monoclonal Antibody
宿主 (Host)	Mouse
应用 (Application)	WB,IHC,IP
种属反应性 (Reactivity)	Human

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	Mouse IgG2b
克隆 (Clonality)	Monoclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	Purified antibody in PBS with 0.05% sodium azide
纯化方式 (Purification)	Affinity Purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	C7orf50
别名 (Alternative Names)	C7orf50
基因 ID (Gene ID)	
蛋白 ID (SwissProt ID)	Q9BRJ6.Purified recombinant fragment of human C7orf50

产品应用 (Application)

稀释比 (Dilution Ratio)	WB 1:500-1:1000,IHC 1:500-1:1000,IP 1:100-1:200
蛋白分子量 (Molecular Weight)	22.1kDa


研究背景 (Background)

产品名称: C7orf50 Mouse Monoclonal Antibody
产品货号: AMM83133




研究领域 (Research Area)


图片 (Image Data)



Western blot analysis using C7orf50 mouse mAb against MCF7 lysate



Lane 1: MCF7 lysate ; Lane 2: C7orf50 immunoprecipitated from MCF7 lysate by AMM83133 ; Lane 3: The same as Lane 2 but use IgG isotype control antibody ;Result: AMM83133 can immunoprecipitate C7orf50;



IHC-P analysis of kidney tissue by C7orf50 antibody (AMM83133)

注意事项 (Note)

产品名称: C7orf50 Mouse Monoclonal Antibody
产品货号: AMM83133



For research use only .