

产品名称: ITGB1BP2 Mouse Monoclonal Antibody  
产品货号: AMM83136



## 产品概述 (Summary)

|                        |                                    |
|------------------------|------------------------------------|
| 产品名称 (Production Name) | ITGB1BP2 Mouse Monoclonal Antibody |
| 描述 (Description)       | Mouse monoclonal Antibody          |
| 宿主 (Host)              | Mouse                              |
| 应用 (Application)       | WB,IP                              |
| 种属反应性 (Reactivity)     | Human, Rat                         |

## 产品性能 (Performance)

|                     |                                                                                          |
|---------------------|------------------------------------------------------------------------------------------|
| 偶联物 (Conjugation)   | Unconjugated                                                                             |
| 修饰 (Modification)   | Unmodified                                                                               |
| 同种型 (Isotype)       | Mouse IgG1                                                                               |
| 克隆 (Clonality)      | Monoclonal                                                                               |
| 形式 (Form)           | Liquid                                                                                   |
| 存放说明 (Storage)      | Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles. |
| 储存溶液 (Buffer)       | Purified antibody in PBS with 0.05% sodium azide                                         |
| 纯化方式 (Purification) | Affinity Purification                                                                    |

## 免疫原信息 (Immunogen)

|                        |                                                        |
|------------------------|--------------------------------------------------------|
| 基因名 (Gene Name)        | ITGB1BP2                                               |
| 别名 (Alternative Names) | ITGB1BP2                                               |
| 基因 ID (Gene ID)        |                                                        |
| 蛋白 ID (SwissProt ID)   | Q9UKP3.Purified recombinant fragment of human ITGB1BP2 |

## 产品应用 (Application)

|                          |                                |
|--------------------------|--------------------------------|
| 稀释比 (Dilution Ratio)     | WB 1:500-1:1000,IP 1:100-1:200 |
| 蛋白分子量 (Molecular Weight) | 38.4kDa                        |


## 研究背景 (Background)

产品名称: ITGB1BP2 Mouse Monoclonal Antibody  
产品货号: AMM83136




## 研究领域 (Research Area)

## 图片 (Image Data)



Western blot analysis using ITGB1BP2 mouse mAb against rat skeletal muscle lysate



Lane 1: Rat skeletal muscle tissue lysate; Lane 2: ITGB1BP2 immunoprecipitated from rat skeletal muscle tissue lysate by AMM83136 ; Lane 3: The same as Lane 2 but use IgG isotype control antibody ;Result: AMM83136 can immunoprecipitate ITGB1BP2;

## 注意事项 (Note)

For research use only .