

产品名称: PRKCB Mouse Monoclonal Antibody
产品货号: AMM83138



产品概述 (Summary)

产品名称 (Production Name)	PRKCB Mouse Monoclonal Antibody
描述 (Description)	Mouse monoclonal Antibody
宿主 (Host)	Mouse
应用 (Application)	WB,IP
种属反应性 (Reactivity)	Human

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	Mouse IgG2b
克隆 (Clonality)	Monoclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	Purified antibody in PBS with 0.05% sodium azide
纯化方式 (Purification)	Affinity Purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	PRKCB
别名 (Alternative Names)	PRKCB
基因 ID (Gene ID)	
蛋白 ID (SwissProt ID)	P05771.Purified recombinant fragment of human PRKCB

产品应用 (Application)

稀释比 (Dilution Ratio)	WB 1:500-1:1000,IP 1:100-1:200
蛋白分子量 (Molecular Weight)	Calculated MW: 77 kDa; Observed MW: 77 kDa


研究背景 (Background)

产品名称: PRKCB Mouse Monoclonal Antibody
产品货号: AMM83138




研究领域 (Research Area)

图片 (Image Data)



Western blot analysis using PRKCB mouse mAb against K-562 lysate



Lane 1: K-562 lysate ; Lane 2: PRKCB immunoprecipitated from K-562 lysate by AMM83138 ; Lane 3: The same as Lane 2 but use IgG isotype control antibody ;Result: AMM83138 can immunoprecipitate PRKCB;

注意事项 (Note)

For research use only .